

Product datasheet for TP307469L

CCDC82 (NM_024725) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human coiled-coil domain containing 82 (CCDC82), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC207469 protein sequence Red =Cloning site Green =Tags(s)


MIHVERRHETRRNSKSHVPEQKSRVDWRRTKRSSISQLLDSDEELDSEEFDSDEELDSDEFENDEELDSN
 KGPDCNKTPGSERELNLSKIQSEGNSKCLINSGNGSTYEETNKKHRNIDLQDQEKHLSQEDNDLNKQ
 TGQIIEDDQEKHLSQEDNDLNKQTGQIIEDDLEEEDIKRGKRRLSSVMCDSDSDSDILVRKVGKRP
 RRVVEDEGSSVMEQKTPEKTLAAQKREKLQKLKELSKQRSRQRSSGRDFEDSEKESCPSSDEVDEEEE
 EDNYEDEDGDDYIIDDFVQDEEGDEENKNQQGEKLTTSQLKLVKRNSLYFSFDHYTHFERVVKALLIN
 ALDESFLGTLTDGTRQKSYAKDMLTSLHYLDNRFVQPRLESLSRSRWKEQYKERVENYSNVSIHLKNPE
 NCSCQACGLHRYCKYSVHLSGELYNTRTMQIDNFMMSHDQVFTVGRICASRTRIYHKLKHFKFELYQECC
 TIAMTEEVEDEQVKETVERIFRRSKENGWIKYGYQLLEEYLNFDYFQEEKFEL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	63.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.


[View online »](#)

RefSeq: NP_079001
 Locus ID: 79780
 UniProt ID: Q8N4S0
 RefSeq Size: 2746
 Cytogenetics: 11q21
 RefSeq ORF: 1632
 Synonyms: HSPC048

Product images:


Coomassie blue staining of purified CCDC82 protein (Cat# [TP307469]). The protein was produced from HEK293T cells transfected with CCDC82 cDNA clone (Cat# [RC207469]) using MegaTran 2.0 (Cat# [TT210002]).