

Product datasheet for TP307481L

C20orf77 (RPRD1B) (NM_021215) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human regulation of nuclear pre-mRNA domain containing 1B (RPRD1B), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC207481 protein sequence Red=Cloning site Green=Tags(s)

MSSFSESALEKKLSELSNSQHSVQTLNLWLIHHRKHAGPIVSVWHRELKAKSNRKLTFLYLANDVIQNS
KRKGPEFTREFESVLVDAFSHVAREADEGCKKPLERLLNIWQERSVYGGFEIQQLKLSMEDSKSPPPKAT
EEKSLKRTFQQIQEEEDDDYPGSYSPQDPSAGPLLTEELIKALQDLENAASGDATVRQKIASLPQEVQD
VSLLEKITDKEAAERLSKTVDEACLLLAEYNGRLAAELEDRRQLARMLVEYTQNQKDVLSKEKKLEEKYK
QKLARVTQVRKELKSHIQSLPDLSSLNVTGGLAPLPSAGDLFSTD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	36.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_067038
Locus ID:	58490



[View online »](#)

UniProt ID:	<u>Q9NQG5</u>
RefSeq Size:	3895
Cytogenetics:	20q11.23
RefSeq ORF:	978
Synonyms:	C20orf77; CREPT; dj1057B20.2; K-H; Kub5-Hera; NET60
Summary:	Interacts with phosphorylated C-terminal heptapeptide repeat domain (CTD) of the largest RNA polymerase II subunit POLR2A, and participates in dephosphorylation of the CTD by RPAP2. Transcriptional regulator which enhances expression of CCND1. Promotes binding of RNA polymerase II to the CCND1 promoter and to the termination region before the poly-A site but decreases its binding after the poly-A site. Prevents RNA polymerase II from reading through the 3' end termination site and may allow it to be recruited back to the promoter through promotion of the formation of a chromatin loop. Also enhances the transcription of a number of other cell cycle-related genes including CDK2, CDK4, CDK6 and cyclin-E but not CDKN1A, CDKN1B or cyclin-A. Promotes cell proliferation.[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified RPRD1B protein (Cat# [TP307481]). The protein was produced from HEK293T cells transfected with RPRD1B cDNA clone (Cat# [RC207481]) using MegaTran 2.0 (Cat# [TT210002]).