

Product datasheet for **TP307482L**

Neuroglobin (NGB) (NM_021257) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human neuroglobin (NGB), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC207482 protein sequence Red =Cloning site Green =Tags(s)
	MERPEPELIRQSWRAVSRSPLEHGTVLFARLFALEPDLLPLFQYNCRQFSSPEDCLSSPEFLDHIRKVML VIDAAVTNVEDLSSLEEYLASLGRKHRAVGVKLSFSTVGESLLYMLEKCLGPAFTPATRAAWSQLYGAV VQAMSRGWDGE
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	16.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_067080
Locus ID:	58157
UniProt ID:	Q9NPG2 , A0M8W9 , F6ZE01 , A0A1L1SVG2
RefSeq Size:	1885




[View online »](#)

Cytogenetics: 14q24.3

RefSeq ORF: 453

Summary: This gene encodes an oxygen-binding protein that is distantly related to members of the globin gene family. It is highly conserved among other vertebrates. It is expressed in the central and peripheral nervous system where it may be involved in increasing oxygen availability and providing protection under hypoxic/ischemic conditions. [provided by RefSeq, Jul 2008]

Product images:



Coomassie blue staining of purified NGB protein (Cat# [TP307482]). The protein was produced from HEK293T cells transfected with NGB cDNA clone (Cat# [RC207482]) using MegaTran 2.0 (Cat# [TT210002]).