

Product datasheet for **TP307497**

GNMT (NM_018960) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins
Description: Recombinant protein of human glycine N-methyltransferase (GNMT), 20 µg
Species: Human
Expression Host: HEK293T
Expression cDNA Clone or AA Sequence: >RC207497 protein sequence
Red=Cloning site **Green**=Tags(s)

MVDSVYRTRSLGVAAEGLPDQYADGEAARVWQLYIGDTRSRTAEYKAWLLGLLRQHGCQRVLDVACGTGV
DSIMLVEEGFSVTSVDASDKMLKYALKERWNRHHEPAFDKWWIEEANWMTLDKDVPQSAEGGFDAVICLG
NSFAHLPDCKGDQSEHRLALKNIASMRAGLLVIDHRNYDHILSTGCAPPGKNIYYKSDLTKDVTTSVL
IVNNKAHMVTLDYTVQVPGAGQDGSPGLSKFRLSYYPHCLASFTELLQAAFGGKCQHSVLGDFKPYKPGQ
TYIPCYFIHVLKRTD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK
Predicted MW: 32.6 kDa
Concentration: >0.05 µg/µL as determined by microplate BCA method
Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage: Store at -80°C.
Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq: [NP_061833](#)
Locus ID: 27232



[View online »](#)

UniProt ID: [Q14749](#)

RefSeq Size: 1091

Cytogenetics: 6p21.1

RefSeq ORF: 885


Synonyms: HEL-S-182mP

Summary: The protein encoded by this gene is an enzyme that catalyzes the conversion of S-adenosyl-L-methionine (along with glycine) to S-adenosyl-L-homocysteine and sarcosine. This protein is found in the cytoplasm and acts as a homotetramer. Defects in this gene are a cause of GNMT deficiency (hypermethioninemia). Alternative splicing results in multiple transcript variants. Naturally occurring readthrough transcription occurs between the upstream CNPY3 (canopy FGF signaling regulator 3) gene and this gene and is represented with GeneID:107080644. [provided by RefSeq, Jan 2016]

Protein Families: Druggable Genome

Protein Pathways: Glycine, serine and threonine metabolism

Product images:



Coomassie blue staining of purified GNMT protein (Cat# TP307497). The protein was produced from HEK293T cells transfected with GNMT cDNA clone (Cat# [RC207497]) using MegaTran 2.0 (Cat# [TT210002]).