

Product datasheet for **TP307537**

DHDH (NM_014475) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human dihydrodiol dehydrogenase (dimeric) (DHDH), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC207537 protein sequence Red =Cloning site Green =Tags(s)

MALRWGIVSVGLISSDFTAVLQTLPRSEHQVVAVAARDLSRAKEFAQKHDIPKAYGSYEELAKDPSVEVA
YIGTQHPQHKAAMVLCCLAAGKAVLCEKPTGVNAAEVREMVAEARSRALFLMEAIWTRFFPASEALRSVLA
QGTLDLRVARAEFGKNLIHVPRAVDRAQAGGALLDIGYCVQFTSMVFGGQKPEKISVVGRRHETGVDD
TVTLLQYPGEVHGSFTCSITVQLSNTASVSGTKGMVQLLNPWCPTLVKGEHKEFPLPPVPKDCNFD
NGAGMSYEAKHVWECLRKGMKESVPIPLSESELLADILEVRKAIGVTFPQDKR

SGPTRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	36.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_055290</u>
Locus ID:	27294



[View online »](#)

UniProt ID: [Q9UQ10](#)

RefSeq Size: 1098

Cytogenetics: 19q13.33


RefSeq ORF: 1002

Synonyms: 2DD; HUM2DD

Summary: This gene encodes an enzyme that belongs to the family of dihydrodiol dehydrogenases, which exist in multiple forms in mammalian tissues and are involved in the metabolism of xenobiotics and sugars. These enzymes catalyze the NADP1-linked oxidation of transdihydrodiols of aromatic hydrocarbons to corresponding catechols. This enzyme is a dimeric dihydrodiol dehydrogenase, and it differs from monomeric dihydrodiol dehydrogenases in its high substrate specificity for trans-dihydrodiols of aromatic hydrocarbons in the oxidative direction. [provided by RefSeq, Jul 2008]

Protein Pathways: Metabolism of xenobiotics by cytochrome P450

Product images:



Coomassie blue staining of purified DHDH protein (Cat# TP307537). The protein was produced from HEK293T cells transfected with DHDH cDNA clone (Cat# [RC207537]) using MegaTran 2.0 (Cat# [TT210002]).