

Product datasheet for **TP307575**

CA7 (NM_001014435) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human carbonic anhydrase VII (CA7), transcript variant 2, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	Recombinant protein was produced with TrueORF clone, RC207575.
Tag:	C-Myc/DDK
Predicted MW:	23.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001014435
Locus ID:	766
UniProt ID:	P43166
RefSeq Size:	1710
Cytogenetics:	16q22.1
RefSeq ORF:	627
Synonyms:	CA-VII; CAVII



[View online »](#)

Summary:


Carbonic anhydrases are a large family of zinc metalloenzymes that catalyze the reversible hydration of carbon dioxide. They participate in a variety of biological processes, including respiration, calcification, acid-base balance, bone resorption, and the formation of aqueous humor, cerebrospinal fluid, saliva, and gastric acid. They show extensive diversity in tissue distribution and in their subcellular localization. The cytosolic protein encoded by this gene is predominantly expressed in the brain and contributes to bicarbonate driven GABAergic neuron excitation. Alternative splicing in the coding region results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Aug 2018]

Protein Families:

Druggable Genome

Protein Pathways:

Nitrogen metabolism

Product images:

Coomassie blue staining of purified CA7 protein (Cat# TP307575). The protein was produced from HEK293T cells transfected with CA7 cDNA clone (Cat# [RC207575]) using MegaTran 2.0 (Cat# [TT210002]).