

Product datasheet for **TP307576M**

APJ Receptor (APLNR) (NM_005161) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human apelin receptor (APLNR), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC207576 protein sequence Red =Cloning site Green =Tags(s)

MEEGGDFDNYYGADNQSECEYTDWKSSGALIPAIYMLVFLGTTGNGLVLWTVFRSSREKRRSADIFIAS
LAVADLTFVTLPLWATYTYRDYDWPFGTFFCKLSSYLIFVNMYSVFCLTGLSFDRYLAIVRPVANARL
RLRVSGAVATAVLWVLAALLAMPVMVLRTTGDLNENTTKVQCYMDYSMVATVSSEWAVEVGLGVSSTTVGF
VVPFTIMLTCYFFIAQTIAGHFRKERIEGLRKRRLLSIIVLVVTFALCWMPYHLVKTLYMLGSLHWP
CDFDLFLMNIFPYCTCISYVNSCLNPFLYAFFDPRFRQACTSMLCCGQSRCAGTSHSSSGEKSASYSSGH
SQGPGPNMGKGGEQMHEKSIPYSQETLVVD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	42.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_005152
Locus ID:	187



[View online »](#)

UniProt ID: [P35414](#), [B3KQN4](#), [B2RDH3](#)

RefSeq Size: 3905

Cytogenetics: 11q12.1

RefSeq ORF: 1140


Synonyms: AGTRL1; APJ; APJR; HG11

Summary: This gene encodes a member of the G protein-coupled receptor gene family. The encoded protein is related to the angiotensin receptor, but is actually an apelin receptor that inhibits adenylate cyclase activity and plays a counter-regulatory role against the pressure action of angiotensin II by exerting hypertensive effect. It functions in the cardiovascular and central nervous systems, in glucose metabolism, in embryonic and tumor angiogenesis and as a human immunodeficiency virus (HIV-1) coreceptor. Two transcript variants resulting from alternative splicing have been identified. [provided by RefSeq, Jul 2009]

Protein Families: Druggable Genome, GPCR, Transmembrane

Protein Pathways: Neuroactive ligand-receptor interaction

Product images:



Coomassie blue staining of purified APLNR protein (Cat# [TP307576]). The protein was produced from HEK293T cells transfected with APLNR cDNA clone (Cat# [RC207576]) using MegaTran 2.0 (Cat# [TT210002]).