

## Product datasheet for **TP307594M**

### ZFP200 (ZNF200) (NM\_003454) Human Recombinant Protein


#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human zinc finger protein 200 (ZNF200), transcript variant 1, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC207594 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	 MMAAKVWPMPPKPKQSFILRVPPDSKLGQDLLRDATNGPKTIHQVLVLEHFLTFLPKPSLVQPSQKVKETL VIMKDVSSSLQNRVHPRPLVKLLPKGVQKEQETVSLYLKANPEELVWFEDLNVFHCQEECVSLDPTQQLT SEKEDDSSVGEMMLLAVNGSNPEGEDPEREPVENEDYREKSSDDDEMDSSLVSQPPDNQEKERLNTSIP QKRKMRNLLVTIENDTPLEELSKYVDISIILTRNRTRRWYTCPLCGKQFNESYLISHQRTHTGKPY DCNHCGKSFNHKTNLNKHRIHTGKPYSCSQCGKNFRQNSHRSRHEGIHIREKIFKCPECCKTFFPKNEE FVLHLQSHEAERPYGCKKCGRRFGRLSNCTRHEKTHSACKTRKQK  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	45.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_003445</a></u>
Locus ID:	7752



[View online »](#)

UniProt ID:	<u>P98182</u> , <u>B3KP91</u>
RefSeq Size:	3387
Cytogenetics:	16p13.3
RefSeq ORF:	1185
Summary:	Could have a role in spermatogenesis.[UniProtKB/Swiss-Prot Function]
Protein Families:	Transcription Factors

**Product images:**

Coomassie blue staining of purified ZNF200 protein (Cat# [TP307594]). The protein was produced from HEK293T cells transfected with ZNF200 cDNA clone (Cat# [RC207594]) using MegaTran 2.0 (Cat# [TT210002]).