

Product datasheet for TP307615

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

GPNMB (NM 001005340) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human glycoprotein (transmembrane) nmb (GPNMB), transcript variant

1, 20 µg

Species: Human Expression Host: HEK293T

Expression cDNA >RC207615 representing NM_001005340

Red=Cloning site Green=Tags(s)

Clone or AA Sequence:

MECLYYFLGFLLLAARLPLDAAKRFHDVLGNERPSAYMREHNQLNGWSSDENDWNEKLYPVWKRGDMRWK NSWKGGRVQAVLTSDSPALVGSNITFAVNLIFPRCQKEDANGNIVYEKNCRNEAGLSADPYVYNWTAWSE DSDGENGTGQSHHNVFPDGKPFPHHPGWRRWNFIYVFHTLGQYFQKLGRCSVRVSVNTANVTLGPQLMEV TVYRRHGRAYVPIAQVKDVYVVTDQIPVFVTMFQKNDRNSSDETFLKDLPIMFDVLIHDPSHFLNYSTIN YKWSFGDNTGLFVSTNHTVNHTYVLNGTFSLNLTVKAAAPGPCPPPPPPPRPSKPTPSLATTLKSYDSNT PGPAGDNPLELSRIPDENCQINRYGHFQATITIVEGILEVNIIQMTDVLMPVPWPESSLIDFVVTCQGSI PTEVCTIISDPTCEITQNTVCSPVDVDEMCLLTVRRTFNGSGTYCVNLTLGDDTSLALTSTLISVPDRDP ASPLRMANSALISVGCLAIFVTVISLLVYKKHKEYNPIENSPGNVVRSKGLSVFLNRAKAVFFPGNQEKD PLLKNQEFKGVS

TRTRPLEQKLISEEDLAANDILDYKDDDDK**V**

Tag: C-Myc/DDK
Predicted MW: 61.5 kDa

Concentration: >0.05 μg/μL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Bioactivity: Pull-down assay (PMID: <u>26751287</u>)

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional

chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience

some loss of protein during the filtration process.





GPNMB (NM_001005340) Human Recombinant Protein - TP307615

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling

conditions. Avoid repeated freeze-thaw cycles.

RefSeq: NP 001005340

Locus ID: 10457
UniProt ID: Q14956
RefSeq Size: 2775
Cytogenetics: 7p15.3
RefSeq ORF: 1716

Synonyms: HGFIN; NMB; PLCA3


Summary: The protein encoded by this gene is a type I transmembrane glycoprotein which shows homology

to the pMEL17 precursor, a melanocyte-specific protein. GPNMB shows expression in the lowly metastatic human melanoma cell lines and xenografts but does not show expression in the highly metastatic cell lines. GPNMB may be involved in growth delay and reduction of metastatic potential. Two transcript variants encoding different isoforms have been found for this gene.

[provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome, Transmembrane

Product images:



Coomassie blue staining of purified GPNMB protein (Cat# TP307615). The protein was produced from HEK293T cells transfected with GPNMB cDNA clone (Cat# [RC207615]) using MegaTran 2.0 (Cat# [TT210002]).