

Product datasheet for **TP307645L**

MRPL47 (NM_177988) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human mitochondrial ribosomal protein L47 (MRPL47), nuclear gene encoding mitochondrial protein, transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	<p>>RC207645 protein sequence Red=Cloning site Green=Tags(s)</p> <p>MAAAGLALLCRRVSSALKSSRSLITPQVPACTGLLHTTLSRKGLEEFFDDPKNWGQEKVKSGAAWTCQQL RNKSNEDLHKLWYVLLKERNMLLTLEQEAKRQRLPMPSPERLDKVDSMDALDKVVQEREDALRLLQTG Q ERARPGAWRRDIFGRIIWHKFKQWVIPWHLNKRYNRKRFFALPYVDHFLRLEREKRARIKARKENLERKK AKILLKKFPHLAEAQKSSLV</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	16.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_817125</u>
Locus ID:	57129


[View online »](#)

UniProt ID:	<u>Q9HD33</u>
RefSeq Size:	1076
Cytogenetics:	3q26.33
RefSeq ORF:	690
Synonyms:	CGI-204; L47mt; MRP-L47; NCM1
Summary:	<p>Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 39S subunit protein. This gene is immediately adjacent to the gene for BAF complex 53 kDa subunit protein a (BAF53a), in a tail-to-tail orientation. Two transcript variants encoding different protein isoforms have been identified. [provided by RefSeq, Jul 2008]</p>

Product images:



Coomassie blue staining of purified MRPL47 protein (Cat# [TP307645]). The protein was produced from HEK293T cells transfected with MRPL47 cDNA clone (Cat# [RC207645]) using MegaTran 2.0 (Cat# [TT210002]).