

## Product datasheet for TP307765L

### PDE4A (NM\_006202) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human phosphodiesterase 4A, cAMP-specific (phosphodiesterase E2 duncce homolog, Drosophila) (PDE4A), transcript variant 4, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC207765 representing NM_006202 Red=Cloning site Green=Tags(s)

MPLVDFFCETCSKPWLVGWWDQFKRMLNRELTHLSEMSRSGNQVSEYISTTFLDKQNEVEIPSPTMKERE  
KQQAPRPRPSQPPPPVPHLQPMSQITGLKMLHMSNSLNSNIPRFGVKTQDELLAQELENLNKWGLNI  
FCVSDYAGGRSLTCIMYMIFQERDLLKKFRIPVDTMVTYMLTLEDHYHADVAYHNSLHAADVQLQSTHVLL  
ATPALDAVFTDLEILAALFAAAIHVDHDPGVSNQFLINTNSELALMYNDESLENHHLAVGFKLLQEDNC  
DIFQNLSKRQRQSLRKMVIDMVLATDMSKHM TLLADLKMVETKKVTSSGVLLLLDNYSDRIQVLRNMVHC  
ADLSNPTKPLELYRQWTD RIMAEFFQQGDRERERGM EISPMCDKHTASVEKSQVGFIDYIVHPLWETWAD  
LVHPDAQEILD TLEDNRDWYSAIRQSPSPPEEESRGP GHPPLPDKFQFELTLEEEEEEEISMAQIPCT  
AQEALTAQGLSGVEEALDATI AWEASPAQESLEVMAQEASLEAELEAVYLTQQAQSTGSAPVAPDEFSSR  
EEFVAVSHSSPSALALQSP LPAWRTLSVSEHAPGLPGLPSTAAEVEAQREHQAAKRAC SACAGTFGED  
TSALPAPGGGGSGGDPT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	72 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Bioactivity:	Binding assay (PMID: <a href="#">27148684</a> )
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.




[View online »](#)

Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_006193</u>
Locus ID:	5141
UniProt ID:	<u>P27815</u>
RefSeq Size:	4255
Cytogenetics:	19p13.2
RefSeq ORF:	1941
Synonyms:	DPDE2; PDE4; PDE46
Summary:	The protein encoded by this gene belongs to the cyclic nucleotide phosphodiesterase (PDE) family, and PDE4 subfamily. This PDE hydrolyzes the second messenger, cAMP, which is a regulator and mediator of a number of cellular responses to extracellular signals. Thus, by regulating the cellular concentration of cAMP, this protein plays a key role in many important physiological processes. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.[provided by RefSeq, Jul 2011]
Protein Families:	Druggable Genome
Protein Pathways:	Progesterone-mediated oocyte maturation, Purine metabolism

### Product images:



PDE4A binds to resveratrol. Purified recombinant PDE4A (2 ug) (OriGene [TP307765]) was incubated with empty (-) or resveratrol-immobilized (+) beads for 4 hours, and the bound PDE4A was detected by Western blotting. The input lane corresponds to recombinant PDE4A (250 ng). Figure cited from Cell Death Dis, PMID: 27148684



Coomassie blue staining of purified PDE4A protein (Cat# [TP307765]). The protein was produced from HEK293T cells transfected with PDE4A cDNA clone (Cat# [RC207765]) using MegaTran 2.0 (Cat# [TT210002]).