

Product datasheet for TP307804

RIBC1 (NM_001031745) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human RIB43A domain with coiled-coils 1 (RIBC1), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC207804 protein sequence Red =Cloning site Green =Tags(s)
	<p>MYNIKQSTDTKEAAAIEARRNREKERQNRFFNVRNRVMGVDVQALNNQVGDRKRREAAERSKEAAYGTSQ VQYDVVVQMLEKEEADRTRQLAKKVQEFREQKQQLKNGREFSLWDPGQVWKGLPTYLSYSNTYPGPASLQ YFSGEDLDRDTRLRMQGGQFRYNLERQQQEQQAKVDENYTDALSNQLRLAMDAQATHLARLEESCAAM MCAMANANKAQAQAVQAGRQRCERQREQKANLAEIQHQSTSDLLTENPQVAQHMPAPYRVLPYCWKGMTPE QQAIRKEQEVQRSKKQHRQAETLDTEWKSQTMSSAQAVLELEEQERELCAVFQRGLGSFNQQLANEQ KAQQDYLSVIYTNQPTAQYHQQFNTSSR</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	43.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001026915</u>



[View online »](#)

Locus ID: 158787
UniProt ID: [Q8N443](#), [A0A024R9X7](#)
RefSeq Size: 1510
Cytogenetics: Xp11.22
RefSeq ORF: 1137
Synonyms: 2610028I09Rik

Product images:



Coomassie blue staining of purified RIBC1 protein (Cat# TP307804). The protein was produced from HEK293T cells transfected with RIBC1 cDNA clone (Cat# [RC207804]) using MegaTran 2.0 (Cat# [TT210002]).