

## Product datasheet for **TP307822M**

### **SYT12 (NM\_177963) Human Recombinant Protein**

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human synaptotagmin XII (SYT12), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC207822 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	<p>MAVDVAEYHLSVIKSPPGWEVGVYAAGALALLGIAAVSLWKLWTSFSPSPFPNYDYRYLQQKYGESC AEAREKRVPAWNAQRASTRGPPSRKGSLSIEDTFESISLGPLLELMGRELDLAPYGTLRKSQSADSLNSI SSVSNTFGQDFTLGQVEVSMEYDTASHTLNVAVMQGKDLLEREEASFESCFMRVSLLPDEQIVGISRIQR NAYSIFDEKFSIPLDPTALEEKSLRFSVFGIDEDERNVSTGVVELKLSVLDLPLQPFSGWLYLQDQNK ADAVGEILLSYLPPTAERLTVVVKAKNLIWTNDKTTADPFVKVYLLQDGRKMSKKKTAVKRDDPNPVF NEAMIFSVPAIVLQDLSLRVTVAESSDGRGDNVGHVIIGPSASGMGTHWNQMLATLRRPVSMWHAVRR N</p> <p><b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b></p>
Tag:	C-Myc/DDK
Predicted MW:	46.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_808878</a></u>




[View online »](#)

Locus ID: 91683  
UniProt ID: [Q8IV01](#), [Q8NDM9](#)  
RefSeq Size: 3541  
Cytogenetics: 11q13.2  
RefSeq ORF: 1263  
Synonyms: SYT11; sytXII

**Summary:** This gene is a member of the synaptotagmin gene family and encodes a protein similar to other family members that mediate calcium-dependent regulation of membrane trafficking in synaptic transmission. Studies of the orthologous gene in rat have shown that the encoded protein selectively modulates spontaneous synaptic-vesicle exocytosis and may also be involved in regulating calcium independent secretion in nonneuronal cells. Alternative splicing results in multiple transcript variants. The gene has previously been referred to as synaptotagmin XI but has been renamed synaptotagmin XII to be standard with mouse and rat official nomenclature.[provided by RefSeq, Apr 2010]

**Protein Families:** Secreted Protein, Transmembrane

### Product images:



Coomassie blue staining of purified SYT12 protein (Cat# [TP307822]). The protein was produced from HEK293T cells transfected with SYT12 cDNA clone (Cat# [RC207822]) using MegaTran 2.0 (Cat# [TT210002]).