

Product datasheet for **TP307852**

CCNJL (NM_024565) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human cyclin J-like (CCNJL), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC207852 protein sequence Red =Cloning site Green =Tags(s)

MMDEPWWEGRVASDVHCTLREKELKLPTFRAHSPLLKSRFFVDILLLSSHCQLCPAARHLAVYLLDHF
MDRYNVTTSKQLYTVAVSCLLLANGVSLSPRLKCSGMISAHCNLHLPGSSNSPASAPHPPTPPQVAET
TGKFEDREDHVPKLEQINSTRILSSQNFTLTKKELLSTELLLEAFSWNLCLPTPAHFLDYLLASVSQK
DHHCHTWPTTCPRKTECLKEYAHYFLEVTLQDHIFYKFQPSVAAACVGASRICLQLSPYWTRDLQRIS
SYSLEHLSTCIEILLWYDNLKDAVAVKSQALAMVPGTPPTPTQVLFQPPAYPALGQPATTLAQFQTPV
QDLCLAYRDSLQAHRSGLLSGSTGSSLHTPYQPLQPLDMCPVVPASLSMHMAIAAEPRHCLATTYGSS
YFSGSHMFPTGCFDR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	48.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_078841</u>



[View online »](#)

Locus ID: 79616
UniProt ID: [Q8IV13](#)
RefSeq Size: 3274
Cytogenetics: 5q33.3
RefSeq ORF: 1305

Product images:



Coomassie blue staining of purified CCNJL protein (Cat# TP307852). The protein was produced from HEK293T cells transfected with CCNJL cDNA clone (Cat# [RC207852]) using MegaTran 2.0 (Cat# [TT210002]).