

Product datasheet for TP307872M

RFX1 (NM_002918) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human regulatory factor X, 1 (influences HLA class II expression) (RFX1), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC207872 protein sequence Red =Cloning site Green =Tags(s)
	<p>MATQAYTELQAAPPSQPPQAPPQAQPQPPPPPPAAPQPPQPPTAAATPQPQYVTELQSPQPQAQPPGG QKQYVTELPVPAPSQPTGAPTSPAPQQYIIVTVSEGAMRASETVSEASPSTASQTGVPTQVWQVQVQ TQQRLLVQTSVQAKPGHVSPQLTNIQVPQQALPTQRLVVQSAAPGSKGGQVSLTVHGTQQVHSPPEQSP VQANSSSSKTAGAPTGTVPQQLQVHGVQSVVPTQERSVQATPQAPKGPVQPLTVQGLQPVHVAQEVQ QLQQVPVPHVYSSQVQVVEGGDASYTASAIRSSTYSYPETPLYTQTASTSYEAAAGTATQVSTPATSQAV ASSGSMMPMYVSGSQVASSASTGAGASNSSGGGGSGGGGGGGGGGGGGSGSTGGGGSGAGTYVIQGGYM LGSASQSYSHTTRASPATVQWLLDNYETAEGVSLPRSTLYCHYLLHCQEQLPEVNAASFGKLIRSVFMG LRTRRLGTRGNSKYHYGLRIKASSPLLRLMEDQQHMAMRGQPFQKQRLKPIQKMEGMTNGVAVGQQPS TGLSDISAQVQYQQLDASRSLPDTFELDLQKVLPEGVGPGDIKAFQVLYREHCEAIVDVMVNLQFTL VETLWKTFRYRNLSPSEAPPLAVHDEAEKRLPKAILVLLSKFEPVLQWTKHCDNVLYQGLVEILIPDVL RPIPSALTQAIRNFAKSLESWLTHAMVNIPEEMLRVKVAAAGAFQAQLRRYTSLNHLAQAARAVLQNTAQ INQMLSDLNRVDFANVQEQASWVCRCEDRVVRLEQDFKVTLQQNSLEQWAAWLDGWSVQLKPYQGSA GFPKAAKLFLLKWSFYSSMVIRDLTLRSAASFGSFHLIRLLYDEYMYLIEHRVAQAKGETPIAVMGEFA NLATSLNPLDPDKDEEEEEEEEESEDELQDISLAAGGESPALGPETLEPPAKLARTDARGLFVQALPSS</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	104.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol



[View online »](#)

Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_002909</u>
Locus ID:	5989
UniProt ID:	<u>P22670</u>
RefSeq Size:	4392
Cytogenetics:	19p13.12
RefSeq ORF:	2937
Synonyms:	EFC; RFX
Summary:	This gene encodes a member of the regulatory factor X (RFX) family of transcription factors, which are characterized by a winged-helix DNA-binding domain. The encoded transcription factor contains an N-terminal activation domain and a C-terminal repression domain, and may activate or repress target gene expression depending on cellular context. This transcription factor has been shown to regulate a wide variety of genes involved in immunity and cancer, including the MHC class II genes and genes that may be involved in cancer progression. This gene exhibits altered expression in glioblastoma and the autoimmune disease systemic lupus erythematosus (SLE). [provided by RefSeq, Jul 2016]
Protein Families:	Druggable Genome

Product images:



Coomassie blue staining of purified RFX1 protein (Cat# [TP307872]). The protein was produced from HEK293T cells transfected with RFX1 cDNA clone (Cat# [RC207872]) using MegaTran 2.0 (Cat# [TT210002]).