

Product datasheet for **TP307890SE**

CD56 (NCAM1) (NM_181351) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human neural cell adhesion molecule 1 (NCAM1), transcript variant 2, secretory expressed in HEK293T cells, 20ug
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC207890 protein sequence Red =Cloning site Green =Tags(s)

MLQTKDLIWTLFFLGTAVSLQVDIVPSQGEISVGESKFFLCQVAGDAKDKDISWFSPNGEKLTPNQQRIS
VWNDDSSSTLIYNANIDDAGIYKCVVTGEDGSESEATVNVKIFQKLMFKNAPTQEFREGEDAVIVCD
VSSLPPTIWKHKGRDVILKKDVRFIVLSNNYLQIRGIKKTDEGTYRCEGRILARGEINFKDIQVIVNV
PPTIQRQNIIVNATANLGQSVTLVCD AEGFPEPTMSWTKDGEQIEQEEDDEKYIFSDSSQLTIKKVDKN
DEAEYICIAENKAGEQDATIHLKVFAPKITYVENQTAMELEEQVTLTCEASGDPIPSITWRTSTRNISS
EEKASWTRPEKQETLDGHMVVRS HARVSSLTKSIQYTDAGEYICTASNTIGQDSQSMYLEVQYAPKLQG
PVAVYTWEGNQVNITCEVFAYPSATISWFRDQQLPSSNYSNIKIYNTPSASYLEVTPDSENDGNYNCT
AVNRIGQESLEFILVQADTPSSPSIDQVEPYSSTAQVQFDEPEATGGVPILKYKAEWRAVGEEVWHSKWY
DAKEASMEGIVTIVGLKPETTYAVRLAALNGKGLGEISAASEFKTQPVGEPSPAPKLEGQMGEDGNSIKV
NLIKQDDGGSPIRHVLRVYRALSSEWKPEIRLPSGSDHVMLKSLDWNAEYEVVVAENQQGKSKAAHFVF
RTSAQPTAIPANGSPTSGLSTGAIVGILIVIFVLLLWVDITCYFLNKCGLFMCIAVNLCGKAGPGAKGK
DMEEGKAAFSKDESKEPIVEVRTEEERTPNHDGGKHTEPNETPLTEPEKGPVEAKPECQETETKPAPAE
VKTVPNDATQTKENENKA


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	96.0 kDa
Concentration:	>50 ug/mL as determined by microplate Bradford method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25mM Tris-HCl, pH7.3, 100mM glycine, 10% glycerol
Note:	For culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.



[View online »](#)

Storage:	Store at -80°C after receiving vials.
Stability:	Stable for at least 1 year from receipt of products under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_851996
Locus ID:	4684
UniProt ID:	P13591
RefSeq Size:	6007
Cytogenetics:	11q23.2
RefSeq ORF:	2574
Synonyms:	CD56; MSK39; NCAM
Summary:	<p>This gene encodes a cell adhesion protein which is a member of the immunoglobulin superfamily. The encoded protein is involved in cell-to-cell interactions as well as cell-matrix interactions during development and differentiation. The encoded protein plays a role in the development of the nervous system by regulating neurogenesis, neurite outgrowth, and cell migration. This protein is also involved in the expansion of T lymphocytes, B lymphocytes and natural killer (NK) cells which play an important role in immune surveillance. This protein plays a role in signal transduction by interacting with fibroblast growth factor receptors, N-cadherin and other components of the extracellular matrix and by triggering signalling cascades involving FYN-focal adhesion kinase (FAK), mitogen-activated protein kinase (MAPK), and phosphatidylinositol 3-kinase (PI3K). One prominent isoform of this gene, cell surface molecule CD56, plays a role in several myeloproliferative disorders such as acute myeloid leukemia and differential expression of this gene is associated with differential disease progression. For example, increased expression of CD56 is correlated with lower survival in acute myeloid leukemia patients whereas increased severity of COVID-19 is correlated with decreased abundance of CD56-expressing NK cells in peripheral blood. Alternative splicing results in multiple transcript variants encoding distinct protein isoforms. [provided by RefSeq, Aug 2020]</p>
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Transmembrane
Protein Pathways:	Cell adhesion molecules (CAMs), Prion diseases

Product images:

Coomassie blue staining of purified NCAM1 protein (Cat #TP307890SE). The protein was produced from mammalian cells transfected with NCAM1 cDNA clone (Cat #[RC207890]).