

Product datasheet for **TP307892L**

PDSS2 (NM_020381) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human prenyl (decaprenyl) diphosphate synthase, subunit 2 (PDSS2), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC207892 protein sequence Red =Cloning site Green =Tags(s)

MNFRQLLLHLPRYLGASGSPRRLWWSPSLDTISSVGSWRGRSSKSPAHWNQWSEAEKIVGYPTSFMSLR
CLLSDELSNIAMQVRKLVGTQHPLLTARGLVHDSWNSLQLRGLVLLISKAAGPSSVNTSCQNYDMVSG
IYSCQRSLAEITELIHIALLVHRGIVNLNELQSSDGPLKDMQFGNKIALLSGDFLLANACNGLALLQNTK
VVELLASALMDLVQGVYHENSTSKESYITDDIGISTWKEQTFLSHGALLAKSCQAAMELAKHDAEVQNMA
FQYGKHMAMSHKINSVQPFIEKETSDSMTFNLNSAPVVLHQEFLGRDLWIKQIREAQEKGRLDYAKLRE
RIKAGKGVTSIDLCRYHGKALEALESFPPSEARSALENIVFAVTRFS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	43.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_065114</u>




[View online »](#)

Locus ID: 57107
UniProt ID: [Q86YH6](#)
RefSeq Size: 3568
Cytogenetics: 6q21
RefSeq ORF: 1197
Synonyms: bA59I9.3; C6orf210; COQ1B; COQ10D3; DLP1; hDLP1

Summary: The protein encoded by this gene is an enzyme that synthesizes the prenyl side-chain of coenzyme Q, or ubiquinone, one of the key elements in the respiratory chain. The gene product catalyzes the formation of all trans-polyprenyl pyrophosphates from isopentyl diphosphate in the assembly of polyisoprenoid side chains, the first step in coenzyme Q biosynthesis. Defects in this gene are a cause of coenzyme Q10 deficiency.[provided by RefSeq, Oct 2009]

Protein Pathways: Terpenoid backbone biosynthesis

Product images:



Coomassie blue staining of purified PDSS2 protein (Cat# [TP307892]). The protein was produced from HEK293T cells transfected with PDSS2 cDNA clone (Cat# [RC207892]) using MegaTran 2.0 (Cat# [TT210002]).