

Product datasheet for TP307905L

alpha 2 Macroglobulin (A2M) (NM_000014) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins
Description: Recombinant protein of human alpha-2-macroglobulin (A2M), 1 mg
Species: Human
Expression Host: HEK293T
Expression cDNA Clone or AA Sequence: >RC207905 protein sequence
Red=Cloning site Green=Tags(s)

MGKMKLLHPSLVLLLLVLLPTDASVSGKPQYMLVPSLLHTETTEKGCVLLSYLNETVTVSASLESVRGN
RSLFTDLEAENDVLHCVAFVAVPKSSNEEVMFLTVQVKGPTQEFKKRTVMVKNEDSLVQVQTKSIYKP
GQTVKFRVSMDFHPLNELIPLVYIQDPKGNRIAQWQSFQLEGGLKQFSFPLSSEPFQGSYKVVVQKK
SGGRTEHPFTVEEFVLPKFEVQVTPKIIITILEEEMNVSVCGLTYGKVPVGHVTVSICRKYSDASDCHG
EDSQAFCEKFSGQLNSHGCFYQQVKTQVFKLKRKEYEMKLHTEAQIQEEGTVELTGRQSSEITRITIKL
SFVKVDShFRQGIPFFGQVRLVDGKGVPIPNKVIFIRGNEANYYSNATTDEHGLVQFSINTTNVMGTSLT
VRVNYKDRSPCYGYQWVSEEHAAHTAYLVFSPSKSFVHLEPMSHELPCGHTQTVQAHYILNGGTLGL
KKLSFYLLIMAKGGIVRTGTHGLLVKQEDMKGHFSISIPVKSIDIAPVARLLIYAVLPTGDVIGDSKYDV
ENCLANKVDLSFSPSQSLPASHAHLRVTAAPQSVCALRAVDQSVLLMKPDAELSASSVYNLLPEKDLTGF
PGPLNDQDDEDCINRHNVIYINGITYTPVSSSTNEKDMYSFLEDMGLKAFTNSKIRKPKMCPQLQYEMHGP
EGLRVGFYESDVMGRGHARLVHVEEPTETVRKYFPETWIWDLVWVNSAGVAEVGVTVPDTITWVWVWV
CLSEDAGLGISSTASLRAFQPFVVELTMPYSVIRGEAFTLKATVLNLYPKCIRVSVQLEASPAFLAVPVE
KEQAPHICANGRQTVSWAVTPKSLGNVNFTVSAEALQSLELQTEVPSVPEHGRKDTVIKPLLVPEGL
EKETTFNSLLCPSGGEVSEELSLKLPNVVEESARASVSVLGDILGSAMQNTQNLQMPYGCGEQNMVLF
APNIYVLDYLNQQLTPEVKSIAIGYLNQYQYRQLNYKHVDGYSYTFGERYGRNQNTWLTAFVLKTF
QARAYIFIDEAHITQALIWLSQRQKDNCGFRSSGSLNNAIKGGVEDEVTLISAYITIALLEIPLTVTHPV
VRNALFCLESAWKTAQEGDHGSHVYTKALLAYAFALAGNQDKRKEVLKSLNNEAVKKNDSVHWERPQPKP
APVGHFYEPQAPSAEVEVMTSYVLLAYLTAQPAPTSEDLSATNIVKWITKQNAQGGFSSTQDTVVALHA
LSKYGAATFTRTGKAAQVTIQSSGTFSSKFQVDNRRLLQVSLPELPGESMKVTEGEGCVLQTSKY
NILPEKEEFPALGVQTLPTQTCDEPKAHTSFQISLSVSYTGSRSASNMAIVDKMVSQFIPLKPTVKMLE
RSNHVSRTEVSSNHVLIYLDKVSNTLSLFFTVLQDVPVRDLKPAIVKVYDYETDEFAIAEYNAPCSKD
LGNA


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK



[View online >](#)

Predicted MW:	160.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_000005
Locus ID:	2
UniProt ID:	P01023
RefSeq Size:	4678
Cytogenetics:	12p13.31
RefSeq ORF:	4422
Synonyms:	A2MD; CPAMD5; FWP007; S863-7
Summary:	The protein encoded by this gene is a protease inhibitor and cytokine transporter. It uses a bait-and-trap mechanism to inhibit a broad spectrum of proteases, including trypsin, thrombin and collagenase. It can also inhibit inflammatory cytokines, and it thus disrupts inflammatory cascades. Mutations in this gene are a cause of alpha-2-macroglobulin deficiency. This gene is implicated in Alzheimer's disease (AD) due to its ability to mediate the clearance and degradation of A-beta, the major component of beta-amyloid deposits. A related pseudogene, which is also located on the p arm of chromosome 12, has been identified. [provided by RefSeq, Nov 2016]
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Complement and coagulation cascades

Product images:

Coomassie blue staining of purified A2M protein (Cat# [TP307905]). The protein was produced from HEK293T cells transfected with A2M cDNA clone (Cat# [RC207905]) using MegaTran 2.0 (Cat# [TT210002]).