

Product datasheet for TP307928L

EHHADH (NM_001966) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human enoyl-Coenzyme A, hydratase/3-hydroxyacyl Coenzyme A dehydrogenase (EHHADH), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC207928 protein sequence Red=Cloning site Green=Tags(s)

MAEYTRLHNALALIRLRNPPVNAISTLLRDIKEGLQKAVIDHTIKAIVICGAEGKFSAGADIRGFSAPR
TFGLTLGHVWDEIQRNEKPVVAIQGMAFGGGLELALGCHYRIAHAEAQVGLPEVTLGLLPGARGTQLLP
RLTGVPAA LDLITSGRRILADEALKLGILDKVNSDPVEEAIRFAQRVSDQPLESRRLCNKPIQSLPNMD
SIFSEALLKMRRQHPGCLAQEACVRAVQAAVQYPYEVGIKKEEELFLYLLQSGQARALQYAFFAERKANK
WSTPSGASWKTASARPVSSVGVVGLGTMGRGIVISFARARIPVIAVSDKNQLATANKMITSVLEKEASK
MQQSGHPWSGPKPRLTSSVKELGGVDLVEAVFEEMSLKKQVFAELSAVCKPEAF LCTNTSALDVDEIAS
STDRPHLVIGTHFFSPA HVMKLEVIP SQYSSPTTIATVMNLSKKIKKIGVVVGNCFGFVGNRMLNPPYIN
QAYFLLEEGSKPEEVDQVLEEF GFKMGPF RVSDLAGLDVGVWKS RKGQGLTGPTLLPGTPARKRGNRRYCP
IPDVLCELGRFGQKTGKGWYQYDKPLGRIHKPDPWLSKFLSRYR KTHHIEPRTISQDEILERCLYSLINE
AFRILGEGIAASPEHIDVVYLHGYGWPRHKGGPMFYASTVGLPTVLEKLQKYRQNP DIPQLEPSDYLLK
LASQGNPPLKEWQSLAGSPSSKL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	79.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.



[View online »](#)

Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001957
Locus ID:	1962
UniProt ID:	Q08426
RefSeq Size:	3870
Cytogenetics:	3q27.2
RefSeq ORF:	2169
Synonyms:	ECHD; FRTS3; L-PBE; LBFP; LBP; MFE1; PBFE
Summary:	The protein encoded by this gene is a bifunctional enzyme and is one of the four enzymes of the peroxisomal beta-oxidation pathway. The N-terminal region of the encoded protein contains enoyl-CoA hydratase activity while the C-terminal region contains 3-hydroxyacyl-CoA dehydrogenase activity. Defects in this gene are a cause of peroxisomal disorders such as Zellweger syndrome. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2009]
Protein Pathways:	beta-Alanine metabolism, Butanoate metabolism, Fatty acid metabolism, Limonene and pinene degradation, Lysine degradation, Metabolic pathways, PPAR signaling pathway, Propanoate metabolism, Tryptophan metabolism, Valine, leucine and isoleucine degradation

Product images:



Coomassie blue staining of purified EHHADH protein (Cat# [TP307928]). The protein was produced from HEK293T cells transfected with EHHADH cDNA clone (Cat# [RC207928]) using MegaTran 2.0 (Cat# [TT210002]).