

Product datasheet for **TP307948**

C12orf60 (NM_175874) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 12 open reading frame 60 (C12orf60), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC207948 protein sequence Red =Cloning site Green =Tags(s)
	 MSSESEKDKERLIQAAKMFFFHVQDLASVINTLTELFSRSMNTQILLMAVRNNSYIKDFFEQMLKIFKEM QSVVDARHDKIQKESLCSKVAMAMCSVVQKSTNVEELHQSAKEVFKSAHTPVIIISVLNSSNILGSLESSL SHLMKFPIMNLQLSDFYTEDTKEQSDVTTSERTRSPPGSSKTTMIDTLKQLDVLKTEDSKNPTKSAADL LEQIVKAMGPILEILQKAIKTMEMNISVFKKASDK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	27.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_787070</u>
Locus ID:	144608
UniProt ID:	<u>Q5U649</u>




[View online »](#)

RefSeq Size: 1139

Cytogenetics: 12p13.1-p12.3

RefSeq ORF: 735

Product images:



Coomassie blue staining of purified C12orf60 protein (Cat# TP307948). The protein was produced from HEK293T cells transfected with C12orf60 cDNA clone (Cat# [RC207948]) using MegaTran 2.0 (Cat# [TT210002]).