

## Product datasheet for **TP308027M**

### **HNRPH3 (HNRNPH3) (NM\_021644) Human Recombinant Protein**

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human heterogeneous nuclear ribonucleoprotein H3 (2H9) (HNRNPH3), transcript variant 2H9A, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC208027 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MDWWMKHNGPNDASDGTVRLRGLPFGCSKEEIVQFFQGLEIVPNGITLTMDYQGRSTGEAFVQFASKEIA ENALGKHKERIGHRYIEIFRSSRSEIKGFYDPPRLLGQRPGPYDRPIGGRGGYAGRGSYGGFDDYGG YNNYGYGNDGFDDMRDGRGMGGHGYGGAGDASSGFHGGHFVHMRGLPFRATENDIANFFSPLNPIRVHI DIGADGRATGEADVFEVTHEDAVAAMSKDKNNMQHRYIELFLNSTPGGGSGMGGSGMGGYGRDGMNQGG YGSVGRMGMGNNYSGGYGTPDGLGGYGRGGGGSGGYGQGGMSGGGWRGMY
	<b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	35.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_067676</a>
Locus ID:	3189



[View online »](#)

UniProt ID: [P31942](#)

RefSeq Size: 2382


Cytogenetics: 10q21.3

RefSeq ORF: 993

Synonyms: 2H9; HNRPH3

**Summary:** This gene belongs to the subfamily of ubiquitously expressed heterogeneous nuclear ribonucleoproteins (hnRNPs). The hnRNPs are RNA binding proteins and they complex with heterogeneous nuclear RNA (hnRNA). These proteins are associated with pre-mRNAs in the nucleus and appear to influence pre-mRNA processing and other aspects of mRNA metabolism and transport. While all of the hnRNPs are present in the nucleus, some seem to shuttle between the nucleus and the cytoplasm. The hnRNP proteins have distinct nucleic acid binding properties. The protein encoded by this gene has two repeats of quasi-RRM domains that bind to RNAs. It is localized in nuclear bodies of the nucleus. This protein is involved in the splicing process and it also participates in early heat shock-induced splicing arrest by transiently leaving the hnRNP complexes. Several alternatively spliced transcript variants have been noted for this gene, however, not all are fully characterized. [provided by RefSeq, Jul 2008]

### Product images:



Coomassie blue staining of purified HNRNPH3 protein (Cat# [TP308027]). The protein was produced from HEK293T cells transfected with HNRNPH3 cDNA clone (Cat# [RC208027]) using MegaTran 2.0 (Cat# [TT210002]).