

Product datasheet for **TP308030M**

BLAP75 (RMI1) (NM_024945) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens RMI1, RecQ mediated genome instability 1, homolog (<i>S. cerevisiae</i>) (RMI1), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA	>RC208030 representing NM_024945
Clone or AA Sequence:	Red=Cloning site Green=Tags(s)

MNVTSIALRAETWLLAAWHVKVPPMWLEACINWIQEENNNVNLSQAQMNKQVFEQWLLTDLRDLEHPLLP
DGILEIPKGELNGFYALQINSLVDVSQPAYSQIQKLRGKNTTNDLVTAEAQVTPKPWEAKPSRMLMLQLT
DGIVQIQGMEYQPIPIHSDLPPGKILYGNISFRLGLVLLKPENVKVLGGEVDALLEEYAEKVLARL
IGEPDLVSVIPNNSNENIPRVTDVLDPALGPSDEELLASLDENDELDTANNDTSSERCFTTGSSSNTIPT
RQSSFEPFVISPRPKEEPSNLSIHVMDGELDDFSLEEALLLEETVQKEQMETKELQPLTFNRNADRSIE
RFSHNPNTTNNFSLTCKNGNNWSEKNVSEQMTNEDKSFSGCPSVRDQNRSIFSVHCNVPLAHDFTNKEKN
LETDNKIKQTSSSDSHSLNNKILNREVNVYQKRNSQISNENDCNLQSCSLRSENSINLSIAMDLYSPP
FVYLSVLMASKPKEVTTVKVKAIVTLTGNLSSSGGIWSITAKVSDGTAYLDVDFVDEILTSLIGFSVPE
MKQSKKDPLQYQKFLEGLQKCQRDLIDLCLMTISFNPSLSKAMVLALQDVNMEHLENLKKRLNK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	70 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_079221
Locus ID:	80010
UniProt ID:	Q9H9A7
RefSeq Size:	3508
Cytogenetics:	9q21.32
RefSeq ORF:	1875
Synonyms:	BLAP75; C9orf76; FAAP75
Summary:	RMI1 is a component of protein complexes that limit DNA crossover formation via the dissolution of double Holliday junctions (Raynard et al., 2006 [PubMed 16595695]).[supplied by OMIM, Mar 2008]

Product images:



Coomassie blue staining of purified RMI1 protein (Cat# [TP308030]). The protein was produced from HEK293T cells transfected with RMI1 cDNA clone (Cat# [RC208030]) using MegaTran 2.0 (Cat# [TT210002]).