

Product datasheet for TP308031M

SMC6L1 (SMC6) (NM_024624) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human structural maintenance of chromosomes 6 (SMC6), transcript variant 2, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC208031 protein sequence Red =Cloning site Green =Tags(s)

MAKRKEENFSSPKNAKRPRQEELEDFDKDGDDECKGTTLTAAEVGIIESIHLKNFMCHSMLGPFKFGSN
VNFVGNNGSGKSAVLTALIVGLGGRVATNRGSSLKGFVKDQNSADISITLRNRGDDAFKASVYGNSI
LIQQHISIDGSRSYKLSATGSVSTRKEELIAILDHFNIQVDNPVSVLTQEMSKQFLQSKNEGDKYKFF
MKATQLEQMKEDESYIMETKERTKEQIHQGEERLTTELKRQCVEKEERFQSIAGLSTMKTNLESLKHEMAW
AVVNEIEKQLNAIRDNIKIGEDRAARLDRKMEEQQVRLNEAEQKYKDIQDKLEKISEETNARAPECMALK
ADVAKKRAYNEAEVLYNRSLNEYKALKKDDQLCKRIEELKKSTDQSLEPERLERQKKISWLKERVKAF
QNQENSVNQEIEQFQQAIEKDKEEHGKIKREELDVKHALSYNQRQLKELKDSKTDRLKRFGNVPALLEA
IDDAYRQGHFTYKPVGPLGACIHLRDPALALAIESCLKGLLQAYCCHNHADERVLQALMKRFYLPGTSRP
PIIWFSEFRNEIYDVRHRAAYHPDFPTVLTALIDNAVANSIDMRGIETVLLIKNNSVARAVMQSQKPP
KNCREAFTADGDQVFAGRYSSENTRPKFLSRDVDSEISDLENEVENKTAQILNLQQHLSALEKDIKHNE
ELLKRCQLHYKELKMKIRKNISEIRELENIIEHQSVDIATLEDEAQENKSKMKMVEEHMEQQKENMEHLK
SLKIEAENKYDAIKFKINQLSELADPLKDELNLADSEVDNQRGKRHYEEKQKEHLDTLNKKKRELDMKE
KELEEKMSQARQICPERIEVEKSASILDKEINRLRQKIQAESHASHGDREEIMRQYQEARETYLDLDSKVR
TLKKFIKLLGEIMEHRFKTYQQFRRCLTLRCKLYFDNLLSQRAYCGKMNFDHKNETLSISVQPGEGNKA
FNDMRALSGGERSFSTVCFILSLWSIAESPFRCLDEFDVMYMDMVNRRIAMDLILKMADSQRFRQFILLTP
QSMSSLPSSKLIRILRMSDPERGQTTLPRPVTQEEDDDQR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	126.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol



[View online »](#)

Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_078900
Locus ID:	79677
UniProt ID:	Q96SB8
RefSeq Size:	5188
Cytogenetics:	2p24.2
RefSeq ORF:	3273
Synonyms:	hSMC6; SMC-6; SMC6L1
Summary:	Core component of the SMC5-SMC6 complex, a complex involved in DNA double-strand breaks by homologous recombination. The complex may promote sister chromatid homologous recombination by recruiting the SMC1-SMC3 cohesin complex to double-strand breaks. The complex is required for telomere maintenance via recombination in ALT (alternative lengthening of telomeres) cell lines and mediates sumoylation of shelterin complex (telosome) components which is proposed to lead to shelterin complex disassembly in ALT-associated PML bodies (APBs). Required for recruitment of telomeres to PML nuclear bodies. SMC5-SMC6 complex may prevent transcription of episomal DNA, such as circular viral DNA genome (PubMed:26983541).[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified SMC6 protein (Cat# [TP308031]). The protein was produced from HEK293T cells transfected with SMC6 cDNA clone (Cat# [RC208031]) using MegaTran 2.0 (Cat# [TT210002]).