

Product datasheet for TP308095

CD53 (NM_000560) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens CD53 molecule (CD53), transcript variant 2, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC208095 representing NM_000560 Red =Cloning site Green =Tags(s) MGMSSLKLLKYVLFNLLFWICGCCILGFGIYLLIHNHNFVLFHNLPSLTLGNVVFVIGSIIMVVAFLG CMGSIKENKCLLMSFFILLIILLAEVTLAILLFVYEQKLNEYVAKGLTDSIHRYHSDNSTKAAWDSIQS FLQCCGINGTSDWTSGPPASCPDRKVEGCYAKARLWFHSNFLYIGIITICVCEVLGMSFALTLCNQI DKTSQTIGL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	24.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_000551</u>
Locus ID:	963



[View online »](#)

UniProt ID: [P19397](#)

RefSeq Size: 1567

Cytogenetics: 1p13.3


RefSeq ORF: 657

Synonyms: MOX44; TSPAN25

Summary: The protein encoded by this gene is a member of the transmembrane 4 superfamily, also known as the tetraspanin family. Most of these members are cell-surface proteins that are characterized by the presence of four hydrophobic domains. The proteins mediate signal transduction events that play a role in the regulation of cell development, activation, growth and motility. This encoded protein is a cell surface glycoprotein that is known to complex with integrins. It contributes to the transduction of CD2-generated signals in T cells and natural killer cells and has been suggested to play a role in growth regulation. Familial deficiency of this gene has been linked to an immunodeficiency associated with recurrent infectious diseases caused by bacteria, fungi and viruses. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2016]

Protein Families: Transmembrane

Product images:



Coomassie blue staining of purified CD53 protein (Cat# TP308095). The protein was produced from HEK293T cells transfected with CD53 cDNA clone (Cat# [RC208095]) using MegaTran 2.0 (Cat# [TT210002]).