

Product datasheet for **TP308117L**

RASL12 (NM_016563) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human RAS-like, family 12 (RASL12), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC208117 protein sequence Red =Cloning site Green =Tags(s)
	MSSVFGKPRAGSGPQSAPLEVNLAAILGRRGAGKSALTVKFLTKRFISEYDPNLEDYSSSEETVDHQPVHL RVMDTADLDTPRNCERYLNWAHAFLVVYSVDSRQSFSSSSYLELLALHAKETQRSIPALLLGKLDMAQ YRQVTKAEGVALAGRFGCLFFEVSACLD FEHVQHVFHEAVREARRELEKSPLTRPLFISEERALPHQAPL TARHGLASCTFNTLSTINLKEMPTVAQAKLVTVKSSRAQSKRKAPTLTLKGFKIF
	TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	29.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_057647</u>
Locus ID:	51285
UniProt ID:	<u>Q9NYN1</u>



[View online »](#)


RefSeq Size: 2603

Cytogenetics: 15q22.31

RefSeq ORF: 798

Synonyms: RIS

Product images:



Coomassie blue staining of purified RASL12 protein (Cat# [TP308117]). The protein was produced from HEK293T cells transfected with RASL12 cDNA clone (Cat# [RC208117]) using MegaTran 2.0 (Cat# [TT210002]).