

Product datasheet for TP308190

FAM78A (NM_033387) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human family with sequence similarity 78, member A (FAM78A), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC208190 protein sequence Red=Cloning site Green=Tags(s)

MGCIQSIGGKARVFREGITVIDVKASIDPVPTSIDESSSVLRYRTPHFRASAQVVMPPIPKKETWVVGW
IQACSHMEFYNYGEGQMSSWELPDLEQEKIQAISDSGVDNYPWYGNNTTCTIVGPTKRDSKFIISMND
NFYPSVTWAVPVSESNVAKLTNIYRDQSFTTWLVATNTSTNDMILQTLHWRMQLSIEVNP NRPLGQRR
LREPIAQDQPKILSKNEPIPPSALVKPNANDAQVLMWRPKYGQPLVIPPKHR


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	31.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_203745
Locus ID:	286336
UniProt ID:	Q5JUQ0 , Q8N2W3



[View online »](#)

RefSeq Size: 3957
Cytogenetics: 9q34.13
RefSeq ORF: 789
Synonyms: C9orf59

Product images:

Coomassie blue staining of purified FAM78A protein (Cat# TP308190). The protein was produced from HEK293T cells transfected with FAM78A cDNA clone (Cat# [RC208190]) using MegaTran 2.0 (Cat# [TT210002]).