

## Product datasheet for TP308296

### C17orf64 (NM\_181707) Human Recombinant Protein


#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 17 open reading frame 64 (C17orf64), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC208296 protein sequence Red=Cloning site Green=Tags(s)
	MLWRFISLFSELEAKQLRRLYKYTKSSQPAKFLVTFCASDAPERSLLADREDSLPKLCHAWGLHSNISGM KERLSNMQTPGQGSPLPGQPRSQDHVKKDSLRELSQKPKLKRKRIKEAPETPETEP
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	14.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_859058</a>
Locus ID:	124773
UniProt ID:	<a href="#">Q86WR6</a>
RefSeq Size:	977
Cytogenetics:	17q23.2


[View online »](#)

RefSeq ORF: 378

### Product images:



Coomassie blue staining of purified C17orf64 protein (Cat# TP308296). The protein was produced from HEK293T cells transfected with C17orf64 cDNA clone (Cat# [RC208296]) using MegaTran 2.0 (Cat# [TT210002]).