

Product datasheet for **TP308300L**

OAF (NM_178507) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human OAF homolog (Drosophila) (OAF), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC208300 protein sequence Red =Cloning site Green =Tags(s)
	MRLPGVPLARPALLLLPLLAPLLGTGAPAE LRVRVRLPDGQVTEESLQADSDADSISLELRKPDGTLVS FTADFKKDVKVFRALILGELEKGQSQFQALCFVTQLQHNEIIPSEAMAKLRQKNPRAVRQAEEVRGLEHL HMDVAVNFSQGALLSPHLHNVC AEVDAIYTRQEDVRFWLEQGVDSSVFEALPKASEQAELPCRQVGD H GKPCVCRYGLSLAWYPCMLKYCHSRDRPTYKCGIRSCQKSYSFDFYVPQRQLCLWDEDPYPG
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	30.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_848602</u>
Locus ID:	220323



[View online »](#)

UniProt ID: Q86UD1
 RefSeq Size: 1957
 Cytogenetics: 11q23.3
 RefSeq ORF: 819
 Synonyms: NS5ATP13TP2

Product images:


Coomassie blue staining of purified OAF protein (Cat# [TP308300]). The protein was produced from HEK293T cells transfected with OAF cDNA clone (Cat# [RC208300]) using MegaTran 2.0 (Cat# [TT210002]).