

Product datasheet for TP308342

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Alpha 1 microglobulin (AMBP) (NM_001633) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human alpha-1-microglobulin/bikunin precursor (AMBP), 20 μg

Species: Human Expression Host: HEK293T

Expression cDNA >RC208342 protein sequence Clone or AA Sequence: Red=Cloning site Green=Tags(s)

MRSLGALLLLLSACLAVSAGPVPTPPDNIQVQENFNISRIYGKWYNLAIGSTCPWLKKIMDRMTVSTLVL GEGATEAEISMTSTRWRKGVCEETSGAYEKTDTDGKFLYHKSKWNITMESYVVHTNYDEYAIFLTKKFSR HHGPTITAKLYGRAPQLRETLLQDFRVVAQGVGIPEDSIFTMADRGECVPGEQEPEPILIPRVRRAVLPQ EEEGSGGGQLVTEVTKKEDSCQLGYSAGPCMGMTSRYFYNGTSMACETFQYGGCMGNGNNFVTEKECLQT CRTVAACNLPIVRGPCRAFIQLWAFDAVKGKCVLFPYGGCQGNGNKFYSEKECREYCGVPGDGDEELLRF

SN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 37.1 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional

chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience

some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling

conditions. Avoid repeated freeze-thaw cycles.

RefSeg: NP 001624

Locus ID: 259



UniProt ID: P02760

RefSeq Size: 1450 Cytogenetics: 9q32 RefSeq ORF: 1056

Synonyms: A1M; EDC1; HCP; HI30; IATIL; ITIL; ITILC; UTI


Summary: This gene encodes a complex glycoprotein secreted in plasma. The precursor is proteolytically

processed into distinct functioning proteins: alpha-1-microglobulin, which belongs to the superfamily of lipocalin transport proteins and may play a role in the regulation of inflammatory processes, and bikunin, which is a urinary trypsin inhibitor belonging to the superfamily of Kunitz-type protease inhibitors and plays an important role in many physiological and pathological processes. This gene is located on chromosome 9 in a cluster of lipocalin genes.

[provided by RefSeq, Jul 2008]

Protein Families: Secreted Protein

Product images:



Coomassie blue staining of purified AMBP protein (Cat# TP308342). The protein was produced from HEK293T cells transfected with AMBP cDNA clone (Cat# [RC208342]) using MegaTran 2.0 (Cat# [TT210002]).