

Product datasheet for **TP308345**

CCDC103 (NM_213607) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human coiled-coil domain containing 103 (CCDC103), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC208345 protein sequence Red =Cloning site Green =Tags(s)
	<p>MERNDIINFKALEKELQAALTADEKYKRENAAKLRAVEQRVASYEEFRGIVLASHLKPLERKDKMGGKRT VPWNCHTIQGRTFQDVATEISPEKAPLQPETSADFYRDWRRHLPSPGPERYQALLQLGGPRLGCLFQTDVVG FGLLGELLVALADHVGPADRAAVLGILCSLASTGRFTLNLSLLSRAERESCKGLFQKLQAMGNPRSVKEG LSWEEQGLEEQSGGLQEEERLLQELLELYQVD</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	27 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_998772
Locus ID:	388389
UniProt ID:	Q8IW40



[View online »](#)

RefSeq Size:	1766
Cytogenetics:	17q21.31
RefSeq ORF:	726
Synonyms:	CILD17; PR46b; SMH
Summary:	This gene encodes a protein that contains a coiled-coil domain. [provided by RefSeq, Apr 2012]

Product images:

Coomassie blue staining of purified CCDC103 protein (Cat# TP308345). The protein was produced from HEK293T cells transfected with CCDC103 cDNA clone (Cat# [RC208345]) using MegaTran 2.0 (Cat# [TT210002]).