

## Product datasheet for TP308352

### RUNDC2B (NM\_001012391) Human Recombinant Protein


#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human RUN domain containing 2B (RUNDC2B), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC208352 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MDEERSSMLPTMAAGPNSILFAINIDNKDLNGQSKFAPTVSDLLKESTQNVTLKESTQGVSSVFREITA SSAISILIKPEQETDPLPVWSRNVSADAKCKKERKKKKQVTNIISFDDEEDEQNSGDMFKKTPGAGESSE DNSDHSSVNIMSAFESPFPGPNSNGSQSSNSWKIDSLNREFGYQKLDVKSIDDEDVDENEDDVGNNSSG RKHRGHSESPEKNGAHSVTQAGVQWHDLGSLQPLPPGFR
	<b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	27 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_001012391</a></u>
Locus ID:	400509
RefSeq Size:	1229



[View online »](#)

Cytogenetics: 16p12.2  
RefSeq ORF: 747  
Synonyms: RUNDC2B; RUNDC2L

**Product images:**

Coomassie blue staining of purified RUNDC2B protein (Cat# TP308352). The protein was produced from HEK293T cells transfected with RUNDC2B cDNA clone (Cat# [RC208352]) using MegaTran 2.0 (Cat# [TT210002]).