

Product datasheet for TP308362M

PTER (NM_001001484) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human phosphotriesterase related (PTER), transcript variant 1, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC208362 protein sequence Red =Cloning site Green =Tags(s) MSSLSGKVQTVLGLVEPSKLGRTLTHEHLAMTFDCCYCPPPPCQEAIKPEIVMKNLYWIKNAYSHKEN LQLNQETEAKEELLYFKANGGGALVENTTTGISRDTQTLKRLAEETGVHISGAGFYVDATHSSETRAM SVEQLTDVLMNEILHGADGTSIKCGIIGEIGCSWPLTESERKVLQATAHAQAQLGCPVIIHPGRSSRAPF QIIRILQEAGADISKTVMShLDRTILDKKELLEFAQLGCYLEYDLFGTELLHYQLGPDIDMPDDNKRIR VRLLEECCEDRILVAHDIHTKTRLMKYGGHGYSHILTNVVPKMLLRGITENVLDKILIENPKQWLTFK TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	38.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001001484</u>
Locus ID:	9317


[View online »](#)

UniProt ID: [Q96BW5](#)
 RefSeq Size: 3826
 Cytogenetics: 10p13
 RefSeq ORF: 1047
 Synonyms: HPHRP; RPR-1

Product images:



Coomassie blue staining of purified PTER protein (Cat# [TP308362]). The protein was produced from HEK293T cells transfected with PTER cDNA clone (Cat# [RC208362]) using MegaTran 2.0 (Cat# [TT210002]).