

Product datasheet for **TP308472L**


PLAC1L (OOSP2) (NM_173801) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human placenta-specific 1-like (PLAC1L), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC208472 protein sequence Red =Cloning site Green =Tags(s)
	MALEVLMLLAVLIWTGAENLHVKISCSLDWLMVSVIPVAESRNLYIFADELHLGMGCPANRIHTYVYEFI YLVRDCGIRTRVVSEETLLFQTELYFTPRNIDHDPQEIHLECSTSRKSVWLTPVSTENEIKLDPSPFIAD FQTAAEELGLLSSPNLL
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	17.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_776162</u>
Locus ID:	219990
UniProt ID:	<u>Q86WS3</u>
RefSeq Size:	1646


[View online »](#)

Cytogenetics: 11q12.1
RefSeq ORF: 474
Synonyms: OOSP2A; PLAC1L; TMEM122
Protein Families: Secreted Protein

Product images:


Coomassie blue staining of purified OOSP2 protein (Cat# [TP308472]). The protein was produced from HEK293T cells transfected with OOSP2 cDNA clone (Cat# [RC208472]) using MegaTran 2.0 (Cat# [TT210002]).