

Product datasheet for **TP308485M**

MOCS2 (NM_004531) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human molybdenum cofactor synthesis 2 (MOCS2), transcript variant 3, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC208485 protein sequence Red =Cloning site Green =Tags(s)
	 MSSLEISSSCFSLETKLPLSPPLVEDSAFEPSPRKDMDEVEEKSKDVINFATAEKLSVDEVSQLVISPLCGA ISLFGTTRNNFEGKKVISLEYEAYLPMAENEVRKICSDIRQKWPVKHIAVFHRLGLVPVSEASIIAVS SAHRAASLEAVSYAIDTLKAKVPIWKKEIYEESSTWKGNKECFWASNS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	20.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_004522
Locus ID:	4338
UniProt ID:	O96007 , A0A024QZS1



[View online »](#)

RefSeq Size: 4200

Cytogenetics: 5q11.2


RefSeq ORF: 564

Synonyms: MCBPE; MOCO1; MOCODB; MPTS

Summary: Eukaryotic molybdoenzymes use a unique molybdenum cofactor (MoCo) consisting of a pterin, termed molybdopterin, and the catalytically active metal molybdenum. MoCo is synthesized from precursor Z by the heterodimeric enzyme molybdopterin synthase. The large and small subunits of molybdopterin synthase are both encoded from this gene by overlapping open reading frames. The proteins were initially thought to be encoded from a bicistronic transcript. They are now thought to be encoded from monocistronic transcripts. Alternatively spliced transcripts have been found for this locus that encode the large and small subunits. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome

Product images:



Coomassie blue staining of purified MOCS2 protein (Cat# [TP308485]). The protein was produced from HEK293T cells transfected with MOCS2 cDNA clone (Cat# [RC208485]) using MegaTran 2.0 (Cat# [TT210002]).