

Product datasheet for **TP308488M**

C6orf62 (NM_030939) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 6 open reading frame 62 (C6orf62), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC208488 protein sequence Red =Cloning site Green =Tags(s)
	 MGDPNSRKKQALNRLRAQLRKKKESLADQFDFKMYIAFVFKEKKKKSALFEVSEVIPVMTNNYEENILKG VRDSSYSLESSLELLQKDWQLHAPRYQSMRRDVGCTQEMDFILWPRNDIEKIVCLLFSRWKESDEPFR PVQAKFEFHGGDYEQFLHVLSRKDKTGIVVNNPNQSVFLFDRQHLQTPKNKATIFKLCSICLYLPQEQ LTHWAVGTIEDHLRPYMPE TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	26.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_112201
Locus ID:	81688
UniProt ID:	Q9GZU0 , A0A024R026



[View online »](#)

RefSeq Size: 2498
Cytogenetics: 6p22.3
RefSeq ORF: 687
Synonyms: dj30M3.2; Nbla00237; XTP12

Product images:

Coomassie blue staining of purified C6orf62 protein (Cat# [TP308488]). The protein was produced from HEK293T cells transfected with C6orf62 cDNA clone (Cat# [RC208488]) using MegaTran 2.0 (Cat# [TT210002]).