

## Product datasheet for **TP308493M**

### SKAP55 (SKAP1) (NM\_001075099) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human src kinase associated phosphoprotein 1 (SKAP1), transcript variant 2, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC208493 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	 MQAAALPEEIRWLLEDAEEFLAEGLRNENLSAVARDHRDHILRGFQQIKARYYWDFQPQGGDIGQDSSDD NHSGTLGLSLTSDAPFLSDYQDEGMEDIVKGAQELDNVIKQGYLEKSKDHSFFGSEWQKRWCVSRGLF YYYANEKSKQPKGTFLIKGYGVRMAPHLRRDSKKESCFELTSQDRRSYEFTATSPAARDWVDQISFLK DLSSLTIPYEEDDEEEKEETYDDIDGFDSPSCGSQCRPTILPGSVGIKEPTEEKEEDIYEVLPDEEHD LEEDESGTRRKGDYASYQGLWDCHGDQPDELSFQRGDLIRILSKEYNMYGWWVGELNSLVGIVPKEYLT TAFEVEER  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	41.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_001068567</a>



[View online »](#)

Locus ID: 8631

UniProt ID: [Q86WV1](#), [V9HW03](#)

RefSeq Size: 1598


Cytogenetics: 17q21.32

RefSeq ORF: 1074

Synonyms: HEL-S-81p; SCAP1; SKAP55

**Summary:** This gene encodes a T cell adaptor protein, a class of intracellular molecules with modular domains capable of recruiting additional proteins but that exhibit no intrinsic enzymatic activity. The encoded protein contains a unique N-terminal region followed by a PH domain and C-terminal SH3 domain. Along with the adhesion and degranulation-promoting adaptor protein, the encoded protein plays a critical role in inside-out signaling by coupling T-cell antigen receptor stimulation to the activation of integrins. [provided by RefSeq, Jul 2008]

### Product images:



Coomassie blue staining of purified SKAP1 protein (Cat# [TP308493]). The protein was produced from HEK293T cells transfected with SKAP1 cDNA clone (Cat# [RC208493]) using MegaTran 2.0 (Cat# [TT210002]).