

Product datasheet for TP308565

FAM83A (NM_032899) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human family with sequence similarity 83, member A (FAM83A), transcript variant 1, 20 µg

Species: Human

Expression Host: HEK293T

Expression cDNA Clone or AA Sequence: >RC208565 protein sequence
Red=Cloning site **Green**=Tags(s)

MSRSRHLGKIRKRLLEDVKSQWVRPARADFSDNESARLATDALLDGGSEAYWRVLSQEGEVDFLSSVEAQY
IQAQAREPPCPPDTLGGAEAGPKGLDSSSLQSGTYFPVASEGSEPALHLSWASAEKPYLKEKSSATVYFQ
TVKHNNIRDLVRRCTRISQVLVILMDVFTDVEIFCDILEAANKRGVFCVLLDQGGVKLFQEMCDKVQI
SDSHLKNISIRSVGEIYCAKSGRKFAGQIREKFIISDWRVFLSGSYSFTWLCGHVHRNLSKFTGQAVE
LFDEEFRHLYASSKPVMGLKSPRLVAPVPPGAAPANGRLSSSSGSASDRTSSNPFSGRSAGSHPGTRSVS
ASSGPCSPAAPHPPPPPRFQPHQGPWGAPSPQAHLSPRPHDGPAAVYSNLGAYRPTRLQLEQLGLVPRL
TPTWRPFLQASPHF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 47.3 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online >](#)

RefSeq: [NP_116288](#)

Locus ID: 84985

UniProt ID: [Q86UY5](#)

RefSeq Size: 4596


Cytogenetics: 8q24.13

RefSeq ORF: 1302

Synonyms: BJ-TSA-9

Summary: Probable proto-oncogene that functions in the epidermal growth factor receptor/EGFR signaling pathway. Activates both RAS/MAPK and PI3K/AKT/TOR signaling cascades downstream of EGFR. Required for the RAS/MAPK signaling cascade activation upon EGFR stimulation, it also activates both signaling cascades independently of EGFR activation. [UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified FAM83A protein (Cat# TP308565). The protein was produced from HEK293T cells transfected with FAM83A cDNA clone (Cat# [RC208565]) using MegaTran 2.0 (Cat# [TT210002]).