

## Product datasheet for **TP308593**

### **TCL1B (NM\_004918) Human Recombinant Protein**

#### **Product data:**


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human T-cell leukemia/lymphoma 1B (TCL1B), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC208593 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MASEASVRLGVPPGRLWIQRPGIYEDEEGRTWVTWVRFNPSRREWARASQGSRYEPSITVHLWQMAVHT RELLSSGQMPFSQLPAVWQLYPRRKYRAADSSFWEIADHGQIDSMEQLVLTYPQPERKD
	<b>TR</b> TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	14.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_004909</a></u>
Locus ID:	9623
UniProt ID:	<u><a href="#">O95988</a></u>
RefSeq Size:	1172



[View online »](#)

Cytogenetics:	14q32.13
RefSeq ORF:	384
Synonyms:	SYN-1; TML1
Summary:	Enhances the phosphorylation and activation of AKT1 and AKT2.[UniProtKB/Swiss-Prot Function]

### Product images:



Coomassie blue staining of purified TCL1B protein (Cat# TP308593). The protein was produced from HEK293T cells transfected with TCL1B cDNA clone (Cat# [RC208593]) using MegaTran 2.0 (Cat# [TT210002]).