

## Product datasheet for **TP308598L**

### **POU6F1 (NM\_002702) Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Recombinant protein of human POU class 6 homeobox 1 (POU6F1), transcript variant 1, 1 mg
<b>Species:</b>	Human
<b>Expression Host:</b>	HEK293T
<b>Expression cDNA Clone or AA Sequence:</b>	>RC208598 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MPGISSQILTNAQGQVIGTLPWVNSASVAAPAPAQSLQVQAVTPQLLLNAQGQVIATLASSPLPPVAV  
RKPSTPESPAKSEVQPIQPTPTVPQPAVVIASPAAPAAKPSASAPIPITCSETPTVSQLVSKPHTPSLDED  
GINLEEIREFAKNFKIRRLSLGLTQTQVGQALTATEGPAYQSAICRFEKLDITPKSAQKLPVLEKWLN  
EAELRNQEGQQNLMFVGGEPSKKRKRRTSFTPQAEALNAYFEKNPLPTGQEITEIAKELNYDREVVVRV  
WFCNRRQTLKNTSKLNVFQIP


**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

<b>Tag:</b>	C-Myc/DDK
<b>Predicted MW:</b>	32.5 kDa
<b>Concentration:</b>	>0.05 µg/µL as determined by microplate BCA method
<b>Purity:</b>	> 80% as determined by SDS-PAGE and Coomassie blue staining
<b>Buffer:</b>	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
<b>Preparation:</b>	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
<b>Note:</b>	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
<b>Storage:</b>	Store at -80°C.
<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<u><a href="#">NP_002693</a></u>
<b>Locus ID:</b>	5463



[View online »](#)

UniProt ID:	<a href="#">Q14863</a>
RefSeq Size:	4715
Cytogenetics:	12q13.13
RefSeq ORF:	903
Synonyms:	BRN5; MPOU; TCFB1
Summary:	Transcription factor that binds preferentially to a variant of the octamer motif (5'-ATGATAAT-3').[UniProtKB/Swiss-Prot Function]
Protein Families:	Transcription Factors

**Product images:**

Coomassie blue staining of purified POU6F1 protein (Cat# [TP308598]). The protein was produced from HEK293T cells transfected with POU6F1 cDNA clone (Cat# [RC208598]) using MegaTran 2.0 (Cat# [TT210002]).