

Product datasheet for **TP308652L**

ZCCHC17 (NM_016505) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human zinc finger, CCHC domain containing 17 (ZCCHC17), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC208652 protein sequence Red =Cloning site Green =Tags(s)
	MNSGRPETMENLPALYTIFQGEVAMVTDYGAFIKIPGCRKQGLVHRTHMSSCRVDKPSEIVDVGDKVWVK LIGREMKNDRIKVSLSMKVVNQGTGKDLDPNNVIIIEQEERRRRSFQDYTGQKITLEAVLNTTCKKCGCKG HFAKDCFMQPGGTKYSLIPDEEEKEEAKSAEFKPDPTRNPSRKRKKEKKKKKHRDRKSSDSDSDSES DTGKRARHTSKDSKAAKKKKKKKKHKKKHKE
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	27.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_057589
Locus ID:	51538
UniProt ID:	Q9NP64 , B4DYA9



[View online »](#)

RefSeq Size:	1635
Cytogenetics:	1p35.2
RefSeq ORF:	723
Synonyms:	HSPC251; pNO40; PS1D

Product images:

Coomassie blue staining of purified ZCCHC17 protein (Cat# [TP308652]). The protein was produced from HEK293T cells transfected with ZCCHC17 cDNA clone (Cat# [RC208652]) using MegaTran 2.0 (Cat# [TT210002]).