

## Product datasheet for **TP308653L**

### CGGBP1 (NM\_003663) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human CGG triplet repeat binding protein 1 (CGGBP1), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC208653 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	 MERFVVTAPPARNRSKTALYVTPDRVTEFGGELHEDGGKLFCTSCNWLNHVRKSAISDHLKSKTHTKR KAEFEEQNVRKKQRPLTASLQCNSTAQTEKVSIVQDFVKMCLEANIPLEKADHPAVRAFLSRHVKNNGSSI PKSDQLRRAYLPDGYENENQLLNSQDC  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	18.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_003654</a>
Locus ID:	8545
UniProt ID:	<a href="#">Q9UFW8</a>



[View online »](#)

RefSeq Size: 4506


Cytogenetics: 3p11.1

RefSeq ORF: 501

Synonyms: CGGBP; p20-CGGBP

**Summary:** This gene encodes a CGG repeat-binding protein that primarily localizes to the nucleus. CGG trinucleotide repeats are implicated in many disorders as they often act as transcription- and translation-regulatory elements, can produce hairpin structures which cause DNA replication errors, and form regions prone to chromosomal breakage. CGG repeats are also targets for CpG methylation. In addition to its ability to bind CGG repeats and regulate transcription, this gene is believed to play a role in DNA damage repair and telomere protection. In vitro studies indicate this protein does not bind to methylated CpG sequences. [provided by RefSeq, Jul 2017]

### Product images:



Coomassie blue staining of purified CGGBP1 protein (Cat# [TP308653]). The protein was produced from HEK293T cells transfected with CGGBP1 cDNA clone (Cat# [RC208653]) using MegaTran 2.0 (Cat# [TT210002]).