

Product datasheet for TP308660M

MTMR14 (NM_001077526) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human myotubularin related protein 14 (MTMR14), transcript variant 1, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC208660 protein sequence Red =Cloning site Green =Tags(s)

MAGARAAAAAASAGSSASSGNQPPQELGLGELLEEFSTQYRAKDGSGTGGSKVERIEKRCLELFGRDYC
FSVIPNTNGDICGHYPRHIVFLEYESSEKEKDTFESTVQVSKLQDLIHRSKMARCRGRFVCPVILFKGKH
ICRSATLAGWGELYGRSGYNYFFSGGADDAWADVEDVTEEDCALRSGDTHLFDKVRGYDIKLLRYLSVKY
ICDLMVENKKVFGMNVTSSEKVDKAQRYADFTLLSIPYPGCEFFKEYKDRDYMAEGLIFNWKQDYVDAP
LSIPDFLTHSLNIDWSQYQCWDLVQQTQNYLKLKLLSLVNSDDDSGLLVHCISGWDRTPLFISLLRSLWA
DGLIHTSLKPTTEILYLTVAIDWFLFGHMLVDRLSKGEEIFFCFNFLKHITSEEFSAKLTQRRKSLPARD
GGFTLEDICMLRRKDRGSTTSLGSDFLVMESSPGATGSFTYEAVELVPAGAPTQAAWRKSHSSSPQSVL
WNRQPSEDRLPSQQGLAEARSSSSSSSNHSDNFFRMGSSPLEVPKPRLAALSDRETRLQEVRSFLAAAY
SSTVGLRAVAPSPSGAIGGLLEQFARGVGLRSISSNAL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	66.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001070994
Locus ID:	64419
UniProt ID:	Q8NCE2 , Q8WYY9
RefSeq Size:	2380
Cytogenetics:	3p25.3
RefSeq ORF:	1794
Synonyms:	C3orf29
Summary:	This gene encodes a myotubularin-related protein. The encoded protein is a phosphoinositide phosphatase that specifically dephosphorylates phosphatidylinositol 3,5-biphosphate and phosphatidylinositol 3-phosphate. Mutations in this gene are correlated with autosomal dominant centronuclear myopathy. Alternate splicing results in multiple transcript variants. A pseudogene of this gene is found on chromosome 18.[provided by RefSeq, Apr 2010]
Protein Families:	Druggable Genome, Phosphatase, Secreted Protein

Product images:



Coomassie blue staining of purified MTMR14 protein (Cat# [TP308660]). The protein was produced from HEK293T cells transfected with MTMR14 cDNA clone (Cat# [RC208660]) using MegaTran 2.0 (Cat# [TT210002]).