

Product datasheet for TP308717M

NUDT10 (NM_153183) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human nudix (nucleoside diphosphate linked moiety X)-type motif 10 (NUDT10), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC208717 protein sequence Red =Cloning site Green =Tags(s) MKCKPNQTRTYDPEGFKKRAACLCFRSEREDEVLVSSSRYPDRWIVPGGGMEPEEEPPGGAAVREYEEA GVKGKLGRLGLGVFEQNQDPEHRTYVYVLTVELLEDWEDSVSIGRKREWFKVEDAIKVLQCHKPVHAEYL EKLKLGGSPTNGNSMAPSSPDSDP TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	18.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_694853</u>
Locus ID:	170685
UniProt ID:	<u>Q8NFP7</u>


[View online »](#)

RefSeq Size:	2018
Cytogenetics:	Xp11.22
RefSeq ORF:	492
Synonyms:	APS2; DIPPP3-alpha; DIPPP3a
Summary:	This gene is a member of the nudix (nucleoside diphosphate linked moiety X)-type motif containing family. The encoded protein is a phosphohydrolase and may regulate the turnover of diphosphoinositol polyphosphates. The turnover of these high-energy diphosphoinositol polyphosphates represents a molecular switching activity with important regulatory consequences. Molecular switching by diphosphoinositol polyphosphates may contribute to the regulation of intracellular trafficking. In some populations putative prostate cancer susceptibility alleles have been identified for this gene. Alternatively spliced transcript variants, which differ only in the 5' UTR, have been found for this gene. [provided by RefSeq, Feb 2015]

Product images:



Coomassie blue staining of purified NUDT10 protein (Cat# [TP308717]). The protein was produced from HEK293T cells transfected with NUDT10 cDNA clone (Cat# [RC208717]) using MegaTran 2.0 (Cat# [TT210002]).