

Product datasheet for **TP308731M**

CTNNA2 (NM_004389) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human catenin (cadherin-associated protein), alpha 2 (CTNNA2), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC208731 protein sequence Red =Cloning site Green =Tags(s)

MTSATSPIILKWDPKSLEIRTLTVERLLEPLVTQVTTLVNTSNKGPSTGKKKGRSKKAHVLAASVEQATQN
FLEKGEQIAKESQDLKEELVAAVEDVRKQGETMRIASSEFADDPCCSSVKRGTMVRAARALLSAVTRLLIL
ADMADVMRLLSHLKIVVEALEAVKNATNEQDLANRFKEFGKEMVKLNYVAARRQELKDPHCRDEMAAAR
GALKKNATMLYTASQAFLRHPDVAATRANRDYVFKQVQEAIAIGISNAAQATSPTDEAKGHTGIGELAAAL
NEFDNKIILDPMTFSEARFRPSLEERLESIIISGAALMADSSCTRDDRERIVAECNAVRQALQDLLSEYM
NNTGRKEKGDPLNIAIDKMTKKTRDLRRQLRKAVMDHISDSFLETNPVLLVIEAAKSGNEKEVKEYAQV
FRESHANKLVEVANLACISNNEEGVKLVRMAATQIDSLCPQVINAALTLAARPOSKVAQDNMDVFKDQWE
KQVRVLTEAVDDITSVDDFLSVSENHILEDVNVKCVIALQEGDVTDLDRTAGAIRGRAARVIHIIAEMEN
YEAGVYTEKVLKLLSETVMRFAEQVEVAIEALSANVPQPFEENEFIDASRLVYDGVDRDIRKAVLMI
RTPEELEDSDFEQEDYDVRSTSQTEDDQLIAGQSARAIMAQLPQEEKAKIAEQVEIFHQEKSKLDAE
VAKWDDSGNDIIVLAKQMCMIMMEMTDFTRGKGPLKNTSDVINAAKKIAEAGSRMDKLARAVADQCPDSA
CKQDLLAYLQRIALYCHQLNICKVKAEVQNLGGELIVSGLDSATSLIQAANKLMNAVVLTKASYVAST
KYQKVYGTAAVNSPVVSWKMKAPEKKPLVKREKPEEFQTRVRRGSQKKHISPVQALSEFKAMDSEF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	100.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



[View online »](#)

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_004380
Locus ID:	1496
UniProt ID:	P26232 , Q49AD3
RefSeq Size:	4005
Cytogenetics:	2p12
RefSeq ORF:	2715
Synonyms:	CAP-R; CAPR; CDCBM9; CT114; CTNR
Summary:	May function as a linker between cadherin adhesion receptors and the cytoskeleton to regulate cell-cell adhesion and differentiation in the nervous system (By similarity). Required for proper regulation of cortical neuronal migration and neurite growth (PubMed:30013181). It acts as negative regulator of Arp2/3 complex activity and Arp2/3-mediated actin polymerization (PubMed:30013181). It thereby suppresses excessive actin branching which would impair neurite growth and stability (PubMed:30013181). Regulates morphological plasticity of synapses and cerebellar and hippocampal lamination during development. Functions in the control of startle modulation (By similarity).[UniProtKB/Swiss-Prot Function]
Protein Families:	Druggable Genome
Protein Pathways:	Adherens junction, Arrhythmogenic right ventricular cardiomyopathy (ARVC), Endometrial cancer, Leukocyte transendothelial migration, Pathways in cancer, Tight junction

Product images:



Coomassie blue staining of purified CTNNA2 protein (Cat# [TP308731]). The protein was produced from HEK293T cells transfected with CTNNA2 cDNA clone (Cat# [RC208731]) using MegaTran 2.0 (Cat# [TT210002]).