

## Product datasheet for **TP308781L**

### LAP2 (TMPO) (NM\_001032283) Human Recombinant Protein

#### Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human thymopoietin (TMPO), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC208781 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	<p> MPEFLEDPSVLTKDKLSELVANNVTLPAGEQRKDVYVQLYLQHLTARNRPPLPAGTNSKGPPDFSSDEE  REPTPVLGSGAAAAGRSRAAVGRKATKKTDKPRQEDKDDLDVTELTNEDLLDQLVKYGVNPGPIVGTTTRK  LYEKKLLKLREQGTESRSSTPLPTISSAENTRQNGSNDSDRYSDNEEDSKIELKLEKREPLKGRAKTPV  TLKQRRVEHNQSYSQAGITETEWTSGSSKGGPLQALTRESTRGSRRTPRKRVTSEHFRIDGPVISESTP  IAETIMASSNESLVNVRTGNFKHASPILPITEFSDIPRRAPKKPLTRAEVGEKTEERRVERDILKEMFP  YEASTPTGISASCRRIKGAAGRPLELSDFRMEESFSSKYVPKYVPLADVKEKTKKGRSIPVWIKILLF  VWVAVFLFLVYQAMETNQVNPFSNFLHVDPRKSN </p> <p><b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b></p>
Tag:	C-Myc/DDK
Predicted MW:	50.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_001027454</a></u>



[View online »](#)

Locus ID:	7112
UniProt ID:	<a href="#">P42167</a>
RefSeq Size:	4186
Cytogenetics:	12q23.1
RefSeq ORF:	1362
Synonyms:	CMD1T; LAP2; LEMD4; PRO0868; TP
Summary:	Through alternative splicing, this gene encodes several distinct LEM domain containing protein isoforms. LEM domain proteins include inner nuclear membrane and intranuclear proteins, and are involved in a variety of cellular functions including gene expression, chromatin organization, and replication and cell cycle control. The encoded alpha isoform is broadly diffuse in the nucleus and contains a lamin binding domain, while the beta and gamma isoforms are localized to the nuclear membrane and contain an HDAC3 interaction domain. The distinct isoforms may compete with each other when acting to chaperone other proteins and regulate transcription. [provided by RefSeq, Aug 2019]
Protein Families:	Stem cell - Pluripotency, Transmembrane

### Product images:



Coomassie blue staining of purified TMPO protein (Cat# [TP308781]). The protein was produced from HEK293T cells transfected with TMPO cDNA clone (Cat# [RC208781]) using MegaTran 2.0 (Cat# [TT210002]).