

Product datasheet for TP308812L

SRGAP1 (NM_020762) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins
Description: Recombinant protein of human SLIT-ROBO Rho GTPase activating protein 1 (SRGAP1), 1 mg
Species: Human
Expression Host: HEK293T
Expression cDNA Clone or AA Sequence: >RC208812 protein sequence
Red=Cloning site Green=Tags(s)

MSTPSRFKKDKKEIIAEYESQVKEIRAQLVEQQKCLEQQTEMRVQLLQDLQDFFRKKAEIETEYSRNLEKL
 AERFMAKTRSTKDHQQYKKDQNLSPVNCWYLLLNQVRRESKDHALTSDIYLNVMRFMQISEDSTRMF
 KKSKEIAFQLHEDLMKVLNELYVMKTYHMYHAESISAESKLKAEKQEEKQIGRSGDPVFHIRLEERHQ
 RRSSVKKIEKMKEKRQAKYSENKLSIKARNEYLLTLEATNASVFKYIHDLSLDLIDCCDLGYHASLNRA
 LRTYLSAEYNLETSRHEGLDIIENAVDNLEPRSDKQRFMEMYPAAFCPPMKFEFQSHMGDEVCQVSAQQP
 VQAEMLRYQQLQSRLATLKIENEVKKTTTEATLQTIQDMVTIEDYDVSECFQHSRSTESVKSTVSETYL
 SKPSIAKRRANQQETEYFYMKLREYLEGSLNITKLQAKHDLLQRTLGEHRAEYMTTRPPNVPPKPKQKH
 RKSRRSQYNTKLFNGDLETFVKDSGQVIPLIVESCIRFINLYGLQHQQGIFRVSGSQVEVNDIKNSFERG
 ENPLADDQSNHDINSVAGVLKLYFRGLENPLFPKERFNDLISCIRIDNLYERALHIRKLLTLPRSVLIV
 MRYLFAFLNHLNQSDENMMDPYNLAICFGPTLMPVPEIQDQVSCQAHVNEIITIIHHETIFPDAKEL
 DGPVYEKCMAGDDYCDSPYSEHGTL EEVDQDAGTEPHTSEDECEPIEAIKFDYVGRSARELSFKKGASL
 LLYHRASEDWWEGRHNGIDGLVPHQYIVVQDMDDTFSDTLSQKADSEASSGPVTEKSSSKDMNSPTDRH
 PDGYLARQRKRGEPPPPVRRPGRGTS DGHCPHPPHALSNSSVDLGSPSLASHPRGLLQNRGLNNDSPERR
 RRPGHGSLTNISRHDSLKKIDSPPIRRSTSSGQYTFNDHKPLDPETIAQDIEETMNTALNELRELERQS
 TAKHAPDVVLDTLEQVKN SPTPATSTESLSPLHNVALRSSEPQIRRSTSSSSDTMSTFKPMVAPRMGVQL
 KPPALRPKPAVLPKTNPTIGPAPPPQGPTDKSCTM

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag: C-Myc/DDK
Predicted MW: 124.1 kDa
Concentration: >0.05 µg/µL as determined by microplate BCA method
Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol



[View online »](#)

Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_065813
Locus ID:	57522
UniProt ID:	Q7Z6B7
RefSeq Size:	8475
Cytogenetics:	12q14.2
RefSeq ORF:	3255
Synonyms:	ARHGAP13; NMTC2
Summary:	The protein encoded by this gene is a GTPase activator, working with the GTPase CDC42 to negatively regulate neuronal migration. The encoded protein interacts with the transmembrane receptor ROBO1 to inactivate CDC42. [provided by RefSeq, Sep 2016]
Protein Pathways:	Axon guidance

Product images:



Coomassie blue staining of purified SRGAP1 protein (Cat# [TP308812]). The protein was produced from HEK293T cells transfected with SRGAP1 cDNA clone (Cat# [RC208812]) using MegaTran 2.0 (Cat# [TT210002]).