

Product datasheet for **TP308857M**

HBQ1 (NM_005331) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human hemoglobin, theta 1 (HBQ1), 100 µg

Species: Human

Expression Host: HEK293T

**Expression cDNA Clone
or AA Sequence:** >RC208857 protein sequence
Red=Cloning site **Green**=Tags(s)

MALSAEDRALVRALWKKLGSNVGVYTTTEALERTFLAFPATKTYFSHLDLSPGSSQVRAHGQKVADALSLA
VERLDDLPHALSALSHLHACQLRVDPAFQLLGHCLLVTLARHYPGDFSPALQASLDKFLSHVISALVSE
YR

TRPLEQKLISEEDLA**ANDILDYKDDDDKV**

Tag: C-Myc/DDK

Predicted MW: 15.3 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_005322](#)

Locus ID: 3049

UniProt ID: [P09105](#), [A0A1K0GUV5](#)

RefSeq Size: 653



[View online »](#)


Cytogenetics: 16p13.3

RefSeq ORF: 426

Synonyms: HBQ

Summary: Theta-globin mRNA is found in human fetal erythroid tissue but not in adult erythroid or other nonerythroid tissue. The theta-1 gene may be expressed very early in embryonic life, perhaps sometime before 5 weeks. Theta-1 is a member of the human alpha-globin gene cluster that involves five functional genes and two pseudogenes. The order of genes is: 5' - zeta - pseudozeta - mu - pseudoalpha-2 - pseudoalpha-1 - alpha-2 - alpha-1 - theta-1 - 3'. Research supports a transcriptionally active role for the gene and a functional role for the peptide in specific cells, possibly those of early erythroid tissue. [provided by RefSeq, Jul 2008]

Product images:



Coomassie blue staining of purified HBQ1 protein (Cat# [TP308857]). The protein was produced from HEK293T cells transfected with HBQ1 cDNA clone (Cat# [RC208857]) using MegaTran 2.0 (Cat# [TT210002]).