

## Product datasheet for TP308867

### ANGPTL3 (NM\_014495) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human angiotensin-like 3 (ANGPTL3), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC208867 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MFTIKLLLFIVPLVISSRIDQDNSSFDLSLSPKSRFAMLDVVKILANGLLQLGHGLKDFVHKTGQIND IFQKLNIFDQSFYDLSLQTSEIKKEEKELRRTTYKLVKNEEVKNMSLELNSKLESLLKILLQKVKY LEEQLTNLIQNQPETPEHPEVTSKTFVEKQDNSIKDLLQTVEDQYKQLNQQHSQIKEIENQLRRTSIQE PTEISLSSKPRAPRTTPFLQLNEIRNVKHDGIPAECTTIYNRGEHTSGMYAIRPSNSQVFHVYCDVISGS PWTLIQHRIDGSQNFNETWENYKYGFRDLGDFWLGLEKIYSIVKQSNYVLRIELEDWVDKNDKHYEYSFY LGNHETNYTLHLVAITGNVPNAIPENKDLVFSTWDHKAKGHFNCPEGYSGGWWWHDECGENNLNGKYNKP RAKSKPERRRGLSWKSQNGRLYSIKSTKMLIHPTDSESEFE
	<b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	51.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_055310</a></u>



[View online »](#)

Locus ID: 27329

UniProt ID: [Q9Y5C1](#)

RefSeq Size: 2951

Cytogenetics: 1p31.3


RefSeq ORF: 1380

Synonyms: ANG-5; ANGPT5; ANL3; FHBL2

**Summary:** This gene encodes a member of a family of secreted proteins that function in angiogenesis. The encoded protein, which is expressed predominantly in the liver, is further processed into an N-terminal coiled-coil domain-containing chain and a C-terminal fibrinogen chain. The N-terminal chain is important for lipid metabolism, while the C-terminal chain may be involved in angiogenesis. Mutations in this gene cause familial hypobetalipoproteinemia type 2. [provided by RefSeq, Aug 2015]

**Protein Families:** Druggable Genome, Secreted Protein

### Product images:



Coomassie blue staining of purified ANGPTL3 protein (Cat# TP308867). The protein was produced from HEK293T cells transfected with ANGPTL3 cDNA clone (Cat# [RC208867]) using MegaTran 2.0 (Cat# [TT210002]).