

Product datasheet for **TP308912M**

HOOK3 (NM_032410) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human hook homolog 3 (Drosophila) (HOOK3), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC208912 protein sequence Red =Cloning site Green =Tags(s)
	MFSVESLERAELCESLLTWIQTFNVDAPCQTVEDLTNGVMAQVLQKIDPAYFDENWLNRIKTEVGDNW R LKISNLKKILKGILDYNHEILGQQINDFTLPDVNLIGEHSDAELGRMLQLILGCAVNCEQKQEYIQA IM MMEESVQHVVMTAIQELMSKESPVSAAGNDAYVDLDRQLKKTTEELNEALSAKEEIAQRCH ELDMQVAAL Q EEKSSLLAENQVLMERLNQSDSIEDPNSPAGRRHLQLQTQLEQLQEETFRLEAAKDDYRIRCEE LEKEIS ELRQQNDELTTLADEAQLKDEIDVLRHSSDKVSKLEGQVESYKKKLEDLGLRRQVKLLEEKNT MYMQN TVSLEELRKANAARSQLETYKRQVVELQNRLESESKADKLD FEYKRLKEKVDLSLQKEKDR LRTERDSL KETIEELRCVQAQEGQLTTQGLMPLGSQESSDSLAAEIVTPEIREKLIRLQHENKMLKLNQEG SDNEKIA LLQSLDDANLRKNELETENRLVNQRLLEVQSQVEELQKSLQDQGSKAEDSVLLKKKLEEHLE KLHEANN ELQKKRAIEDLEPRFNNSSLKIEELQEALRKKEEEMKQMEERYKKYLEKAKSVIRLDPKQ NQGAPEI QALKNQLQERDRLFHSLEKEYEKTQSQREMEEKYIVSAWYNMGMTLHKKAAEDRLASTGSGQ SFLARQR Q ATSSRRSYPGHVQPATAR TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	82.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



[View online »](#)

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_115786
Locus ID:	84376
UniProt ID:	Q86VS8
RefSeq Size:	14460
Cytogenetics:	8p11.21
RefSeq ORF:	2154
Synonyms:	HK3
Summary:	Hook proteins are cytosolic coiled-coil proteins that contain conserved N-terminal domains, which attach to microtubules, and more divergent C-terminal domains, which mediate binding to organelles. The Drosophila Hook protein is a component of the endocytic compartment.[supplied by OMIM, Apr 2004]

Product images:



Coomassie blue staining of purified HOOK3 protein (Cat# [TP308912]). The protein was produced from HEK293T cells transfected with HOOK3 cDNA clone (Cat# [RC208912]) using MegaTran 2.0 (Cat# [TT210002]).