

Product datasheet for **TP308954**

PXDC1 (NM_183373) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 6 open reading frame 145 (C6orf145), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC208954 protein sequence Red =Cloning site Green =Tags(s)
	 MASAVFEGTSLVNMFVRGCVWNGIRRLIVSRRGDEEEFFEIRTEWSDRSVLYLHRSLADLGRWLWQRLRDA FPEDRSELAQGPLRQGLVAIKEAHDIELRLNEVEKLLKTIISMPCKYSRSEVLTFFERSPLDQVLKNDN VHKIQPSFQSPVKISEIMRSNGFCLANTETIVIDHSIPNGRDQHLGVDPTHELFENGSEFPSELEDGDDP AAYVTNLSYYHLVPFETDIWD TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	26.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_899229
Locus ID:	221749
UniProt ID:	Q5TGL8



[View online »](#)

RefSeq Size: 2211
Cytogenetics: 6p25.2
RefSeq ORF: 693
Synonyms: C6orf145

Product images:

Coomassie blue staining of purified PXDC1 protein (Cat# TP308954). The protein was produced from HEK293T cells transfected with PXDC1 cDNA clone (Cat# [RC208954]) using MegaTran 2.0 (Cat# [TT210002]).